

Article code	7981845 (M90800345)
Article name acc. DoC	Freka Tube FR 10 linear, 80 cm, ENFit
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11553	Solvent	●CHN (Cyclohexanone)
C22501	Aceton - for analysis	●Acetone (2-Propanone, Dimethyl ketone)
C22572	Lubricant	●PDMS (Polydimethylsiloxane)
C22704	Fixation tape	●polyester fleece non-woven ●polyacrylate adhesive ●siliconized paper liner
C30728	Tube Tip FR 10	●PUR (Polyurethane) with 10% BaSO ₄ ●Lubricant (Amidwachs) ●Slip additive Colorslip 103 100 % carboxylic acid ester
C30749	ENFit male 10 FR	●PET (Polyethyleneterephthalate) (Copolyester) ●colourant violet
C30823	ENFit female cap	●HDPE (high density Polyethylene) ●colourant violet
C30824	Strap for ENFit cap, ring-Ø 10.5 mm	●LDPE (low density polyethylene), transparent
C30906	Tube 2.4 x 3.4 mm	●PUR (Polyurethane) + PUR with 60% BaSO ₄ + Colourant (white)
C31526	Thinner for printing ink	●1-Methoxy-2-Propanol (2-Methoxymethylethoxy) propanol
C31527	Printing ink, black	●Stamping ink
C31534	Tube printing ink	●Butanone ●Solvent Black 29 ●Ethanol ●Cellulose ●Isopropyle Acetate ●2-Naphtole
C31535	Tube printing ink thinner	●Butanone ●Acetone
Packaging		
	Label sticker 145 x 94 mm	●80 g/sqm adhesive paper (adhesive based on acrylic)
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 249 x 148 x 100 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
	Automat reel without print	●Paper
	Paper 148 x 280 mm	●Paper

Article code	Article name	Material
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide