

|                                |                                                              |
|--------------------------------|--------------------------------------------------------------|
| <b>Article code</b>            | 8081271                                                      |
| <b>Article name acc. DoC</b>   | Ambix Intrastick® Safe 20Gx17mm max 21 bar/300 psi-max 4ml/s |
| <b>Class acc. MDR 2017/745</b> | Ila                                                          |

| Article code | Article name                  | Material                                                                                             |
|--------------|-------------------------------|------------------------------------------------------------------------------------------------------|
| C11515       | Solvent                       | ●MEK (Methylethylketone)                                                                             |
| C11553       | Solvent                       | ●CHN (Cyclohexanone)                                                                                 |
| C22044       | Handle Piece IS 20G           | ● MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)                                           |
| C22046       | TELESKOPIC SEGMENT 1          | ● PC (Polycarbonate) based on Bisphenole A                                                           |
| C22047       | TELESKOPIC SEGMENT 2          | ●PC (Polycarbonate)                                                                                  |
| C22048       | TELESKOPIC SEGMENT 3          | ●PC (Polycarbonate) based on Bisphenole A                                                            |
| C22558       | Glue                          | ●ECA (Ethyl-2-Cyanoacrylate) Hydrochinon (0.01- < 0.1% w/w)                                          |
| C23139       | Cannula 20G X 29,5 MM         | ●Stainless Steel 1.4301(chrome-nickel steel)                                                         |
| C23148       | Adhesive Plate, printed 20G   | ●polyester fleece<br>●polyacrylate adhesive<br>●siliconized paper liner                              |
| C30762       | Protection Cap                | ●LDPE (low density polyethylen)                                                                      |
| C30969       | Tube 3.5 x 0.58 mm            | ●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO4 + Colourant (white)                        |
| C30993       | Tube                          | ●PUR (Polyurethane)                                                                                  |
| C31771       | Bottom Plate                  | ●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)                                            |
| J1040008     |                               |                                                                                                      |
| C22141       | holding plate                 | ●PP (Polypropylene) contains SEBS (Styrene-ethylene/ butylene-styrene)                               |
| C90476       | Cawiton PR7723C-11 Z2         | ●(SEBS) Styrol ethylene butylene styrol compound with carrier PP (polypropylene) Bormed RD 804 CF-11 |
| J1040245     |                               |                                                                                                      |
| C31536       | Silicone fluid MED-420-350 cP | ●Liquid copolymer                                                                                    |
| C31636       | K-Zero Female Luer            | ●Eastman Tritan Copolyester MX731                                                                    |
| C31637       | K-Zero Male Luer              | ●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)<br>●colour batch light blue, 3%           |
| C31639       | K-Zero Microbial Barrier      | ●Liquid Silicone Rubber<br>●Phenyl Oil                                                               |
| Packaging    |                               |                                                                                                      |

|                      |                                     |                                                                                               |
|----------------------|-------------------------------------|-----------------------------------------------------------------------------------------------|
|                      | Label sticker, size: 104 x 148 mm   | ●Paper for laserprinting                                                                      |
|                      | Label sticker, size: 97 x 67.7 mm   | ●Paper 70 g / qm<br>●adhesive (Acrylate-dispersion)                                           |
|                      | IFU Ambix Intrastick                | ●Paper                                                                                        |
|                      | Paper Tyvek blank                   | ●HDPE (high density polyethylene)                                                             |
|                      | Self-adhesive tape, white           | ●PP (Polypropylene)<br>●adhesive: natural rubber                                              |
|                      | Natron Paper 1000 x 750 mm          | ●Color brown                                                                                  |
|                      | Stretch film                        | ●LLDPE (linear low density polyethylene)                                                      |
|                      | Stretch film foil plate             | ●LDPE (low density polyethylene)                                                              |
|                      | Cardboard 586 x 251 x 165 mm        | ●Corrugated cardboard, color brown<br>●Printing FK Logo: Blue- Pantone 300U                   |
|                      | LDPE- Bag Tube 400 X 60 mm          | ●LDPE (low density polyethylen)                                                               |
|                      | Protective Layer 200 x 370 mm       | ●cardboard                                                                                    |
|                      | Folding box 245 x 87 x 46 mm        | ●Corrugated cardboard, color brown inside and outside                                         |
|                      | Unprinted Tyvek paper               | ●Paper                                                                                        |
|                      | PETG-Foil                           | ●PETG (Polyethyleneterephthalate Glycol)-Foil                                                 |
|                      | Black ink                           | ●Printing ink                                                                                 |
|                      | Orange ink                          | ●Printing ink                                                                                 |
| <b>Sterilization</b> |                                     |                                                                                               |
|                      | Indicator tape for EO sterilization | ●Process indicator tape for Ethylene Oxide (EO) sterilization                                 |
|                      | Sterilization gas, Poland           | ●100% Ethylene oxide, liquefied grade one (GAZ 4-134)<br>●alternative: > 99.9% Ethylene oxide |