

|                                |                                                               |
|--------------------------------|---------------------------------------------------------------|
| <b>Article code</b>            | 8084221                                                       |
| <b>Article name acc. DoC</b>   | Ambix NONCORstick® Safe 19Gx17mm max 21 bar/300 psi-max 5ml/s |
| <b>Class acc. MDR 2017/745</b> | Ila                                                           |

| Article code     | Article name                      | Material                                                                                                                      |
|------------------|-----------------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| C11515           | Solvent                           | ●MEK (Methylethylketone)                                                                                                      |
| C11553           | Solvent                           | ●CHN (Cyclohexanone)                                                                                                          |
| C12579           | Identification Label 19G          | ●Liner: PE (Polyethylene) film with transparent adhesive<br>●Coating: UV Varnish<br>●Printing Ink: HP ElectroInk Black Q4133B |
| C21941           | Tube clamp with centering, white  | ●POM (Polyoxymethylene)<br>●Color Carrier: resin-ester polymer with pigment dispersion                                        |
| C22043           | Handle Piece IS 19G               | ●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)                                                                     |
| C22046           | TELESKOPIC SEGMENT 1              | ●PC (Polycarbonate)                                                                                                           |
| C22047           | TELESKOPIC SEGMENT 2              | ●PC (Polycarbonate)                                                                                                           |
| C22048           | TELESKOPIC SEGMENT 3              | ●PC (Polycarbonate)                                                                                                           |
| C22558           | Glue                              | ●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon                                                                                     |
| C23134           | Cannula 19G X 29,5 MM             | ●Stainless Steel                                                                                                              |
| C30730           | Holding Segment AMBIX NONCOR      | ●PP (Polypropylene)                                                                                                           |
| C30732           | Skin protection layer PAD, AMBIX  | ●polyester fleece<br>●polyacrylate adhesive<br>●siliconized paper liner                                                       |
| C30761           | Protection Cap                    | ●LDPE (low density polyethylen)                                                                                               |
| C30969           | Tube 3.5 x 0.58 mm                | ●PUR (Polyurethane) and PUR (Polyurethane) with BaSO4 + Colourant (white)                                                     |
| C30993           | Tube                              | ●PUR (Polyurethane)                                                                                                           |
| C31641           | Bottom Plate                      | ●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)                                                                     |
| <b>J1040245</b>  |                                   |                                                                                                                               |
| C31536           | Silicone fluid MED-420-350 cP     | ●Liquid copolymer                                                                                                             |
| C31636           | K-Zero Female Luer                | ●Eastman Tritan Copolyester MX731                                                                                             |
| C31637           | K-Zero Male Luer                  | ●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)<br>●colour batch light blue                                        |
| C31639           | K-Zero Microbial Barrier          | ●Liquid Silicone Rubber<br>●Phenyl Oil                                                                                        |
| <b>Packaging</b> |                                   |                                                                                                                               |
|                  | Label sticker, size: 104 x 148 mm | ●Paper for laserprinting                                                                                                      |
|                  | Label sticker, size: 97 x 67.7 mm | ●Paper 70 g / qm<br>●adhesive (Acrylate-dispersion)                                                                           |
|                  | IFU AMBIX NONCORSTICK SAFE        | ●Printed IFU                                                                                                                  |
|                  | Self-adhesive tape, white         | ●PP (Polypropylene)<br>●adhesive: natural rubber                                                                              |
|                  | Natron Paper 1000 x 750 mm        | ●Color brown                                                                                                                  |

| Article code         | Article name                        | Material                                                                            |
|----------------------|-------------------------------------|-------------------------------------------------------------------------------------|
|                      | Stretch film                        | ●LLDPE (linear low density polyethylene)                                            |
|                      | Stretch film foil plate             | ●LDPE (low density polyethylene)                                                    |
|                      | Cardboard 586 x 251 x 165 mm        | ●Corrugated cardboard, color brown<br>●Printing FK Logo: Blue- Pantone 300U         |
|                      | Protective Layer 200 x 370 mm       | ●cardboard                                                                          |
|                      | Cardboard 249 x 148 x 100 mm        | ●Corrugated cardboard, color brown<br>●Printing FK Logo: Blue- Pantone 300U         |
|                      | PETG-Foil                           | ●PETG (Polyethyleneterephthalate Glycol)-Foil                                       |
|                      | Corrugated Cardboard 245MM          | ●Corrugated Cardboard                                                               |
|                      | TYVEK Paper NONCORSTICK             | ●Printed Label (raw material: Tyvek white pantone 021U<br>variable data: WLK660070) |
| <b>Sterilization</b> |                                     |                                                                                     |
|                      | Indicator tape for EO sterilization | ●Process indicator tape for Ethylene Oxide (EO)<br>sterilization                    |
|                      | Sterilization gas                   | ●Ethylene oxide                                                                     |