

Article code	8084241
Article name acc. DoC	Ambix NONCORstick® Safe 19Gx27mm max 21 bar/300 psi-max 5ml/s
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C12579	Identification Label 19G	●Liner: PE (Polyethylene) film with transparent adhesive ●Coating: UV Varnish ●Printing Ink: HP ElectroInk Black Q4133B
C21941	Tube clamp with centering, white	●POM (Polyoxymethylene) ●Color Carrier: resin-ester polymer with pigment dispersion
C22043	Handle Piece IS 19G	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C22046	TELESKOPIC SEGMENT 1	●PC (Polycarbonate)
C22047	TELESKOPIC SEGMENT 2	●PC (Polycarbonate)
C22048	TELESKOPIC SEGMENT 3	●PC (Polycarbonate)
C22558	Glue	●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon
C23136	Cannula 19G X 39,5 MM	●Stainless Steel
C30730	Holding Segment AMBIX NONCOR	●PP (Polypropylene)
C30732	Skin protection layer PAD, AMBIX	●polyester fleece ●polyacrylate adhesive ●siliconized paper liner
C30761	Protection Cap	●LDPE (low density polyethylen)
C30969	Tube 3.5 x 0.58 mm	●PUR (Polyurethane) and PUR (Polyurethane) with BaSO4 + Colourant (white)
C30993	Tube	●PUR (Polyurethane)
C31641	Bottom Plate	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
J1040245		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31637	K-Zero Male Luer	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●colour batch light blue
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	Label sticker, size: 97 x 67.7 mm	●Paper 70 g / qm ●adhesive (Acrylate-dispersion)
	IFU AMBIX NONCORSTICK SAFE	●Printed IFU
	Self-adhesive tape, white	●PP (Polypropylene) ●adhesive: natural rubber
	Natron Paper 1000 x 750 mm	●Color brown

Article code	Article name	Material
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 586 x 251 x 165 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Protective Layer 200 x 370 mm	●cardboard
	Cardboard 249 x 148 x 100 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	PETG-Foil	●PETG (Polyethyleneterephthalate Glycol)-Foil
	Corrugated Cardboard 245MM	●Corrugated Cardboard
	TYVEK Paper NONCORSTICK	●Printed Label (raw material: Tyvek white pantone 021U variable data: WLK660070)
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide