

Article code	8084251
Article name acc. DoC	Ambix NONCORstick® Safe 19Gx35mm max 21 bar/300 psi-max 5ml/s
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C12579	Indentification Label 19G	●Liner: PE (Polyethylene) film with transparent adhesive ●Coating: UV Varnish ●Printing Ink: HP ElectroInk Black Q4133B
C21941	Tube clamp with centering, white	●POM (Polyoxymethylene) ●Color Carrier: resin-ester polymer with pigment dispersion
C22046	TELESKOPIC SEGMENT 1	●PC (Polycarbonate) based on Bisphenole A
C22047	TELESKOPIC SEGMENT 2	●PC (Polycarbonate)
C22048	TELESKOPIC SEGMENT 3	●PC (Polycarbonate) based on Bisphenole A
C22558	Glue	●ECA (Ethyl-2-Cyanoacrylate) + Hydrochinon
C30732	Skin protection layer PAD, AMBIX NO	●polyester fleece ●polyacrylate adhesive ●siliconized paper liner
C30848	TELESCOPIC SEGMENT 4	●PC (Polycarbonate) based on Bisphenole A
C30849	Holding Segment	● PP (Polypropylene)
C30851	Handle Piece 19G	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C30854	Cannula 19G X 47,5 mm	●Stainless Steel 1.4301
C30858	Bottom Plate	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C30969	Tube 3.5 x 0.58 mm	●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO4 + Colourant (white)
C30993	Tube	●PUR (Polyurethane)
C31527	Printing ink, black	●Stamping ink
C31734	Needle protection cap 1,4 x 2,3 x 46 mm (creme)	● LDPE
J1040245		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31637	K-Zero Male Luer	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) ●colour batch light blue

Article code	Article name	Material
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		
C11600	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
C11608	Label sticker, size: 97 x 67.7 mm	●Paper 70 g / qm ●adhesive (Acrylate-dispersion)
C12820	IFU Ambix NONCORstick	●Paper
C22621	Self-adhesive tape, white	●PP (Polypropylene) ●adhesive: natural rubber
C22645	Natron Paper 1000 x 750 mm	●Color brown
C22687	Reels for automatic packaging unprinted Tyvek paper	●N/A
C31045	Stretch film	●LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	●LDPE (low density polyethylene)
C31061A	Cardboard 586 x 251 x 165 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
C31077	Protective Layer 200 x 370 mm	●cardboard
C31092A	Cardboard 249 x 148 x 100 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
C31351	Compound foil 30/70 µm	●PA (Polyamide)/PE (Polyethylene)
C31526	Thinner for printing ink	●1-Methoxy-2-Propanol (2-Metholymethylethoxy) propanol
C64868	Unprinted Tyvek paper	●Tyvek Paper
Sterilization		
C11525	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
C31500	Sterilization gas	●Ethylene oxide