

Article code	8084271
Article name acc. DoC	Ambix NONCORstick® Safe 20Gx17mm max 21 bar/300 psi-max 4ml/s
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C12578	Identification Label 20G	●Liner: PE(Polyethylene) film with transparent adhesive ●Coating: UV Varnish ●Printing Ink: HP ElectroInk Black
C22044	Handle Piece IS 20G	● MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C22046	TELESKOPIC SEGMENT 1	● PC (Polycarbonate) based on Bisphenole A
C22047	TELESKOPIC SEGMENT 2	●PC (Polycarbonate)
C22048	TELESKOPIC SEGMENT 3	●PC (Polycarbonate) based on Bisphenole A
C22558	Glue	●ECA (Ethyl-2-Cyanoacrylate) Hydrochinon (0.01- < 0.1% w/w)
C23139	Cannula 20G X 29,5 MM	●Stainless Steel 1.4301(Chrome-nickel steel)
C30730	Holding Segment AMBIX NONCOR	●PP(Polypropylene)
C30732	Skin protection layer PAD, Ambix NO	●polyester fleece ●polyacrylate adhesive ●siliconized paper liner
C30733	Tube Clamp, GE	●PP(Polypropylene)
C30762	Protection Cap	●LDPE (low density polyethylen)
C30969	Tube 3.5 x 0.58 mm	●PUR (Polyurethane) and PUR (Polyurethane) with 60% BaSO4 + Colourant (white)
C30993	Tube	●PUR (Polyurethane)
C31771	Bottom Plate	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
J1040245		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31637	K-Zero Male Luer	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil

Packaging		
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	Label sticker, size: 97 x 67.7 mm	●Paper 70 g / qm ●adhesive (Acrylate-dispersion)
	IFU Ambix NONCORstick	●Paper
	Paper Tyvek blank	●HDPE (high density polyethylene)
	Self-adhesive tape, white	●PP (Polypropylene) ●adhesive: natural rubber
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 586 x 251 x 165 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	Protective Layer 200 x 370 mm	●cardboard
	Cardboard 249 x 148 x 100mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue-Pantone 300U
	Unprinted Tyvek paper	●Paper
	PETG-Foil	●PETG (Polyethyleneterephthalate Glycol)-Foil
	Black ink	●Printing ink
	Orange ink	●Printing ink
	Corrugated Cardboard 245 x 145mm	●3-layer cardboard
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	●100% Ethylene oxide, liquefied grade one (GAZ 4-134) ●alternative: > 99.9% Ethylene oxide