

Article code	8088041
Article name acc. DoC	Therastick 28G 12 mm 60 cm
Class acc. Directive 93/42/EEC	Ila

Article code	Article name	Material
C10542	Protective cap	●HDPE (high density polyethylene)
C10669	LL neg adaptor, moulded to tube	●EVA (ethylene vinyl acetate copolymer)
C10953	Tube 0.42 x 0.5 mm	●LDPE (low density polyethylene)
C11145N	Plaster 39 x 44,8 mm	●single-coated nonwoven ●Base release paper ●Silicone coating
C11511	Lubricant	●Aminofunctional Siloxane
C11512	Solvent	●Naphtha (petroleum), hydrotreated light
C21633	Cannula 28G, 0.36 mm OD x 28.5 mm, stich lenght 12 ± 1 mm	●Stainless steel
C21714	Handle piece	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C21718	Needle holder , moulded to tube	●LDPE (low density polyethylene)
Packaging		
	Carton 400 x 220 x 380 mm	●Corrugated cardboard, Flute B
	Folding box 210 x 190 x 87 mm	●Cardboard
	Adhesive tape paper	●gummed paper tapes Polyeste filament yarn for strengthening
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	IFU	●Paper
	Label sticker 105 x 144 mm	●Paper
	Paper Tyvek blank	●HDPE (high density polyethylene)
	Foil 400 µm	●A-PET (amorphous -Polyethyleneterephthalate)
	PE banderole 38 x 10 x 0.12 mm	●LDPE (low density polyethylene)
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Tyvek paper	●printed Tyvek paper
Sterilization		
	Sterilization gas	●Ethylene oxide
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization