

Article code	8502822
Article name acc. DoC	Closing-Cap Luer, yellow
Class acc. Directive 93/42/EEC	Is

Article code	Article name	Material
C10616	Closing Cap, yellow	●ABS (Acrylonitrile Butadiene Styrene) ●Carrier: ABS (Acrylonitrile Butadiene Styrene) - Components: Inorganic and organic pigments
Packaging		
	Carton 400 x 220 x 380 mm	●Corrugated cardboard, Flute B
	Carton box 185 x 100 x 72 mm	●Cardboard
	Adhesive tape paper	●gummed paper tapes - Polyeste filament yarn for strengthening
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	Label sticker, size: 97 x 67.7 mm	●Paper 70 g / qm ●adhesive (Acrylate-dispersion)
	White paper for IFU printing	●100% wood pulp
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Peel paper	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide