

Article code	9000782LDT
Article name acc. DoC	Freka® Connect ENFit/ProNeo LDT 2.5 ml
Class acc. Directive 93/42/EEC	Is

Article code	Article name	Material
Material 153	Syringe barrel	•PP (Polypropylene)
Material 154	Syringe plunger	•PP (Polypropylene) •coulorant (violet) RAL 4005
Material 155	Syringe gasket	•IR (Isoprene Rubber)
Material 156	Syringe lubricant	•Silicone oil
Packaging		
Material 157	Print	•PPE Ink Series
Material 160	Secondary packaging (30 syringes in their primary packaging)	•Paper
Material 161	Tertiary packaging (20 boxes of 30 syringes)	•corrugated paper
	Primary packaging	•Porous Material (medical grade glue paper) •PA/PE film •print ink (water based)
Sterilization		
Material 158	Sterilisation	•ETO (Ethylene oxide)