

Article code	9000783LDT
Article name acc. DoC	Freka® Connect 5 ml ENFit/ProNeo LDT (30x1)
Class acc. Directive 93/42/EEC	Ila

Article code	Article name	Material
Material 153	Syringe barrel	●PP (Polypropylene)
Material 154	Syringe plunger	●PP (Polypropylene) ●colorant (violet) RAL 4005
Material 155	Syringe gasket	●IR (Isoprene Rubber)
Material 156	Syringe lubricant	●Silicone oil
Packaging		
Material 157	Print	●PPE Ink Series
Material 160	Secondary packaging (30 syringes in their primary packaging)	●Paper
Material 162	Tertiary packaging (16 boxes of 30 syringes)	●corrugated paper
	Primary packaging	●Porous Material (medical grade glue paper) ●PA/PE film ●print ink (water based)
Sterilization		
Material 158	Sterilisation	●ETO (Ethylene oxide)