

Article code	9000785S
Article name acc. DoC	Freka® Connect 20ml ENFit/ProNeo (30x1)
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
	Syringe Barrel	Polypropylene
	Syringe Plunger	Polypropylene
	Syringe Gasket	Isoprene Rubber
	Syringe Lubricant	Silicone Oil
	Graduation Print	PPE Ink Series
Packaging		
	Primary Container	Porous Material (Medical grade glue paper) with Polypropylene & Polyethylene (Medical composite packaging film)
	Secondary Container	Hard Paper
	Tertiary Container	Corrugated Paper
	Pallet	Euro Pallet 1200 x 800mm
	Wrapping	Shrink
Sterilization		
	Sterilization gas	ETO (Ethylene Oxide)