

Article code	9004122
Article name acc. DoC	PE 100 cm
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10535	Protective cap lock	●LDPE (low density polyethylene)
C10542	Protective cap	●HDPE (high density polyethylene)
C10945	Tube 1.0 x 2.0 mm	●LDPE (low density polyethylene)
C30086	Male connector LL	●HDPE (high density polyethylene)
C30093	Female connector LL	●HDPE (high density polyethylene)
Packaging		
C11020	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
C11032	HDPE-Seitenfaltenbeutel	●HDPE (high density polyethylene)
C11054A	Carton box	●Cardboard
C11119	Flat pouch 145 x 255 x 0.08 mm	●LDPE (low density polyethylene) / LLDPE (linear low density polyethylene)
C11301	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
C11517	Adhesive tape paper	●gummed paper tapes Polyeste filament yarn for strengthening
C11600	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
C11945	White paper for IFU printing	●100% wood pulp
C31045	Stretch film	●LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	●LDPE (low density polyethylene)
C61365	Paper 66.7 x 107.5 mm	●Paper
Sterilization		
C11525	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
C31500	Sterilization gas	●Ethylene oxide