

Article code	9004222
Article name acc. DoC	PVC, 200 cm, light-protected, yellow
Class acc. MDR 2017/745	IIa

Article code	Article name	Material
C10511	protective cap LL	● PS (polystyrene) with butadien HIPS (high impact polystyrene) ● PS (polystyrene) with silicone ; ● PS (Polystyrene)
C10636	Female connector LL	● PVC (Polyvinylchloride) without plasticizer
C10914	Tube 1.5 x 2.7 mm	● PVC (Polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate), UV-absorber
C11553	Solvent	● CHN (Cyclohexanone)
J1010900		
C10505	Protective cap lock	● PP (Polypropylene)
C10603	male connector Luer Lock	● MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene) ● MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene),
Packaging		
	Tubular foil for packaging	● Film 500m x 350 mm x 40 µm ● Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Carton box	● Cardboard
	Compound foil 70 µm	● PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Adhesive tape paper	● gummed paper tapes Polyeste filament yarn for strengthening
	PE banderole 45 x 25 x 0.05 mm	● HDPE (high density polyethylene) banderole
	White paper for IFU printing	● 100% wood pulp
	Roll Label 105 x 148 mm	● Paper (blank white) and Adhesive
	Stretch film	● LLDPE (linear low density polyethylene)
	Stretch film foil plate	● LDPE (low density polyethylene)
	Paper 66.7 x 150 mm	● Paper
Sterilization		
	Indicator tape for EO sterilization	● Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	● 100% Ethylene oxide, liquefied grade one (GAZ 4-134) ● alternative: > 99.9% Ethylene oxide