

Article code	9004442
Article name acc. DoC	Injectomat Line PVC 250 cm
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11553	Solvent	CHN (Cyclohexanone),
C30937S	PVC-Tube 1,5x2,7	PVC (Polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
J1010900		
C10505	Protective cap lock	PP (Polypropylene)
C10603	male connector Luer Lock	MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene) MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)
J1010906		
C10511	protective cap LL	PS (polystyrene) with butadien HIPS (high impact polystyrene) (polystyrene) with silicone PS (Polystyrene),
C10636	Female connector LL	PVC (Polyvinylchloride) without plasticizer
Packaging		
	Carton 400x220x380 mm	carton
	Carton 400 x 220 x 380 mm	Corrugated cardboard, Flute B
	Tubular foil for packaging	Film 500m x 350 mm x 40 µm Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Flat pouch 145 x 255 x 0.08 mm	LDPE (low density polyethylene) / LLDPE (linear low density polyethylene)
	Compound foil 70 µm	PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Adhesive tape paper	gummed paper tapes Polyester filament yarn for strengthening
	PE banderole 45 x 25 x 0.05 mm	HDPE (high density polyethylene) banderole
	Label sticker, size: 104 x 148 mm	Paper for laserprinting
	White paper for IFU printing	100% wood pulp
	Roll Label 105 x 148 mm	Paper (blank white) and Adhesive
	PE bag with bottomseal, transparent	LDPE (low density polyethylene) Film thickness: 45 µm
	Stretch film	LLDPE (linear low density polyethylene)
	Stretch film foil plate	LDPE (low density polyethylene)
	PE (Polyethylene) film broadband	PE (Polyethylene), white
	Paper 66.7 x 150 mm	Paper

Article code	Article name	Material
	Paper Injectomat Line PVC 250 cm	Paper
Sterilization		
	Sterilization gas, Nanchang	Ethylene oxide (90%) + Carbon dioxide (10%)
	Indicator tape for EO sterilization	Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	100% Ethylene oxide, liquefied grade one (GAZ 4-134) alternative: > 99.9% Ethylene oxide,