

|                                       |                            |
|---------------------------------------|----------------------------|
| <b>Article code</b>                   | 9014342                    |
| <b>Article name acc. DoC</b>          | Injectomat Line PVC 200 cm |
| <b>Class acc. Directive 93/42/EEC</b> | IIa                        |

| <b>Article code</b> | <b>Article name</b>                 | <b>Material</b>                                                                                                                     |
|---------------------|-------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| C11553              | Solvent                             | CHN (Cyclohexanone),                                                                                                                |
| C30937S             | PVC-Tube 1,5x2,7                    | PVC (Polyvinylchloride) with plasticizer<br>TEHTM (Tri-(2-ethylhexyl) trimellitate)                                                 |
| <b>J1010900</b>     |                                     |                                                                                                                                     |
| C10505              | Protective cap lock                 | PP (Polypropylene)                                                                                                                  |
| C10603              | male connector Luer Lock            | MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)<br>MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)                |
| <b>J1010906</b>     |                                     |                                                                                                                                     |
| C10511              | protective cap LL                   | PS (polystyrene) with butadien HIPS (high impact polystyrene)<br>(polystyrene) with silicone<br>PS (Polystyrene),                   |
| C10636              | Female connector LL                 | PVC (Polyvinylchloride) without plasticizer                                                                                         |
| <b>Packaging</b>    |                                     |                                                                                                                                     |
|                     | Carton 400x220x380 mm               | carton                                                                                                                              |
|                     | Carton 400 x 220 x 380 mm           | Corrugated cardboard, Flute B                                                                                                       |
|                     | Tubular foil for packaging          | Film 500m x 350 mm x 40 µm<br>Mixture of different Polyolefines PE (Polyethylene),<br>LDPE (Polyethylene, low density), transparent |
|                     | Flat pouch 145 x 255 x 0.08 mm      | LDPE (low density polyethylene) / LLDPE (linear low density polyethylene)                                                           |
|                     | Compound foil 70 µm                 | PA (Polyamide) / adhesive layer / PE (Polyethylene)                                                                                 |
|                     | Adhesive tape paper                 | gummed paper tapes Polyester filament yarn for strengthening                                                                        |
|                     | PE banderole 45 x 25 x 0.05 mm      | HDPE (high density polyethylene) banderole                                                                                          |
|                     | Label sticker, size: 104 x 148 mm   | Paper for laserprinting                                                                                                             |
|                     | White paper for IFU printing        | 100% wood pulp                                                                                                                      |
|                     |                                     |                                                                                                                                     |
|                     | PE bag with bottomseal, transparent | LDPE (low density polyethylene)<br>Film thickness: 45 µm                                                                            |
|                     | Stretch film                        | LLDPE (linear low density polyethylene)                                                                                             |
|                     | Stretch film foil plate             | LDPE (low density polyethylene)                                                                                                     |
|                     | PE (Polyethylene) film broadband    | PE (Polyethylene), white                                                                                                            |
|                     | PE band, Length: 500 +/-            | ●HDPE (high density polyethylene)                                                                                                   |

| Article code         | Article name                                | Material                                                                                     |
|----------------------|---------------------------------------------|----------------------------------------------------------------------------------------------|
|                      | 5m,<br>Width: 45 +/- 3 m, 50 µm<br>+/-10 µm |                                                                                              |
|                      | Paper                                       | Paper                                                                                        |
|                      | Paper "Injectomat Line                      | Paper                                                                                        |
| <b>Sterilization</b> |                                             |                                                                                              |
|                      | Sterilization gas,<br>Nanchang              | Ethylene oxide (90%) + Carbon dioxide (10%)                                                  |
|                      | Indicator tape for EO<br>sterilization      | Process indicator tape for Ethylene Oxide (EO)<br>sterilization                              |
|                      | Sterilization gas, Poland                   | 100% Ethylene oxide, liquefied grade one (GAZ 4-134)<br>alternative: > 99.9% Ethylene oxide, |