

Material List
Finished Product: F2859271



Internal code	Designation	Polymer
330288	Protective cap, male, white	●ABS (Acrylonitrile Butadiene Styrene)
J1040332		
C11553	Solvent	●CHN (Cyclohexanone), CAS-no: 108-94-1
C31740	Luer Lock connector female, white	●ABS (Acrylonitrile Butadiene Styrene)
C31747	Clamp, white	●PP (Polypropylene)
C31749	4-way Y-Connector	● PVC (Polyvinylchloride), plasticizer TEHTM (Tri-2-(Ethylhexyl) Trimellitate)
C31754	Spike Assembly, Filter 3µm, Cap	●Spike: ABS (Acrylonitrile Butadiene Styrene); ●Filter 3 µm: PE (Polyethylene) / Versapor 3000RC; ●Cap: LDPE (Low-Density Polyethylene);
C60963	Tube 4.8 x 6.8	●PVC (Polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
<u>Packaging</u>		
330048	Label sticker 105 x 74 mm	●Paper, white
340010	Sterile bag 220 x 300 mm	●PET (Polyethyleneterephthalate) ●PP (Polypropylene) ●Paper
341020	Carton OD 400 x 400 x 250 mm	●corrugated cardboard
342014	Label sticker 105 x 148 mm	●Paper
350-04	Transfer film	●Transfer film, black
<u>Sterilization</u>		
322036	Sterilization indicator, brown	●Sterilization indicator ETO
Material 240	Sterilization	●ETO (Ethylene Oxide)