

Article code	F2860291
Article name acc. DoC	Freka® Mix Bag 2000 mL +6-filling leads (Twist Off)
Class acc. MDR 2017/745	Is

Article code	Article name	Material
1	Tube foil 0.35 mm	EVA (ethylene vinyl acetate copolymer)
2	Coex tube 6.3 x 7.95 x 40 mm	EVA (ethylene vinyl acetate copolymer) / PUR (Polyurethane)
3	Coex tube 4.9 x 6.5 mm	EVA (ethylene vinyl acetate copolymer) / PUR (Polyurethane)
4	Breakable injection site with protective cap	ABS (Acrylonitrile Butadiene Styrene) PP (Polypropylene) Polyisoprene (Latex-free)
5	Non-reopening clamp, nature	PP (Polypropylene)
6	Locking cap female, blue	ABS (Acrylonitrile-Butadiene-Styrene) Colourant
7	LL connector male, white	ABS (Acrylonitrile Butadiene Styrene)
8	Handle board, blue	PP (Polypropylene), Colourant
9	Injection port for Breakable Injection Site + Protective cap for injection port	PC (Polycarbonate)/PP (Polypropylene), Polyisoprene rubber LDPE (Low-density Polyethylene)
10	Twist-off membrane	PVC (Polyvinylchloride), DEHP-free
11	Bag imprint	Stamp foil blue
12	Solvent	THF (Tetrahydrofuran)
13	Solvent	CHN (Cyclohexanone)
J1040333		
C11553	Solvent	CHN (Cyclohexanone)
C31740	Luer Lock connector female, white	ABS (Acrylonitrile Butadiene Styrene)
C31747	Clamp, white	PP (Polypropylene)
C31749	4-way Y-Connector	PVC (Polyvinylchloride), plasticizer TEHTM (Tri-2-(Ethylhexyl) Trimellitate)
C31750	3-way Y-Connector	PVC (Polyvinylchloride)
C31754	Spike Assembly, Filter 3µm, Cap	Spike: ABS (Acrylonitrile Butadiene Styrene); Filter 3 µm: PE (Polyethylene) / Versapor 3000RC; Cap: LDPE (Low-Density Polyethylene);
C60963	Tube 4.8 x 6.8	PVC (Polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
J1031759	Assembled Drip Chamber	Infusion drip chamber: PVC (Polyvinylchloride); 2-way cover: PVC (Polyvinylchloride); 25 µm disk filter: PA (Polyamide); Blue latex-free injection site: PP (Polypropylene) / PI

Article code	Article name	Material
		(Polyisoprene) / ABS (Acrylonitrile-Butadiene-Styrene)
Packaging		
	Sterile gas paper	Paper
	Compound foil 80 µm	PA (Polyamide) / PE (Polyethylene)
	Foil band 180 x 100 mm	HDPE (high density polyethylene)
	Carton OD 400 x 400 x 350 mm	corrugated cardboard
	Label sticker 105 x 148 mm	Paper
	Instructions for use Freka Mix	Paper
	Transfer film	Transfer film, black
	Transfer print	Transfer foil, black
	Flexoprint	Printing ink, black
Sterilization		
	Sterilization indicator, brown	Sterilization indicator ETO
	Sterilization	ETO (Ethylene Oxide)