

IHE Integration Statement

Vendor	Product Name	Version	Date
Fresenius Kabi	Vigilant Software Suite	1.4	Nov. 2022

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Consistent Time	Time Client	None
Device Enterprise Communication	Device Observation Reporter	None
Infusion Pump Event Communication	Device Observation Reporter	None

Internet address for vendor's IHE information

https://key2.fresenius-kabi.com/sites/default/files/document_public/

Links to Standards Conformance Statements for the Implementation

HL7	Not Applicable
DICOM	Not Applicable

Link to general information on IHE

www.ihe.net