

Article code	M20830020
Article name acc. DoC	Ambix Intrastick® Safe 22G x 24mm max 21 bar/300 psi-max 1,5ml/s
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C22045	Handle Piece IS 22G	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C22046	TELESKOPIC SEGMENT 1	●PC (Polycarbonate)
C22047	TELESKOPIC SEGMENT 2	●PC (Polycarbonate)
C22048	TELESKOPIC SEGMENT 3	●PC (Polycarbonate)
C22558	Glue	●ECA (Ethyl-2-Cyanoacrylate); Hydrochinon
C23149	Adhesive Plate, printed 22G	●polyester fleece ●polyacrylate adhesive ●siliconized paper liner
C30761	Protection Cap	●LDPE (low density polyethylen)
C30969	Tube 3.5 x 0.58 mm	●PUR (Polyurethane) and PUR (Polyurethane) with BaSO4 + Colourant (white)
C30993	Tube	●PUR (Polyurethane)
C31641	Bottom Plate	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C31738	Port needle 22G 36,5 mm	●Stainless steel
J1040008		
C22141	holding plate	●PP (Polypropylene) contains SEBS (Styrene- ethylene/ butylene-styrene)
C90476	Cawiton PR7723C-11 Z2	●(SEBS) Styrol ethylene butylene styrol compound with carrier PP (polypropylene)
J1040245		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31637	K-Zero Male Luer	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●colour batch light blue
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	Label sticker, size: 97 x 67.7 mm	●Paper 70 g / qm ●adhesive (Acrylate-dispersion)
	IFU Ambix Intrastick	●Paper
	Paper Tyvek blank	●HDPE (high density polyethylene)
	Self-adhesive tape, white	●PP (Polypropylene) ●adhesive: natural rubber
	Natron Paper 1000 x 750 mm	●Color brown
	Stretch film	●LLDPE (linear low density polyethylene)

Article code	Article name	Material
	Stretch film foil plate	●LDPE (low density polyethylene)
	Cardboard 586 x 251 x 165 mm	●Corrugated cardboard, color brown ●Printing FK Logo: Blue- Pantone 300U
	LDPE-Bag Tube 400 X 60 mm	●LDPE (low density polyethylen)
	Protective Layer 200 x 370 mm	●cardboard
	Folding box 245 x 87 x 46 mm	●Corrugated cardboard, color brown inside and outside
	Unprinted Tyvek paper	●Paper
	PETG-Foil	●PETG (Polyethyleneterephthalate Glycol)-Foil
	Black ink	●Printing ink
	Orange ink	● Printing ink
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide