

Article code	M46421910
Article name acc. DoC	Ambix nova stationary Set
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10898	Filter Speed Flow 1,2 µm	● Housing: MBS (Methylacrylate / butadiene/styrene) ● fluid filter: 1.2 µm hydrophilic PES Polyethersulfone) ● air filter: 0.03 µm hydrophobic PTFE (Polytetrafluoroethylene)
C11503	Solvent	● THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11553	Solvent	● CHN (Cyclohexanone)
C30936	Tube 3.0 x 4.1 mm	● PVC (Polyvinylchloride) soft, noDOP, 34% TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30949	Tube 2.5 x 4.22 mm	● PVC (Polyvinylchloride)
J1010110		
C10178	Drip chamber upper part 21 E	● PS (polystyrene), ● Colourant white, 50 % TiO₂ in HIPS (high impact polystyrene)
C10221	moulded ring white for drip chamber 21	● (98%) PS (Polystyrene), ● (2%) Colourant, carrier material PS (Polystyrene), contains 50% TiO₂
C10265	Drip chamber lower part 21 V92	● SBS (Styrene Butadiene Copolymer)
C10345	air ventilation flap	● (50%) CAS-No.: 9003-07-0 1-propene homopolymer: PP (Polypropylene) ● (50%) CAS-No.: 9010-79-1 1-propene copolymer (heterophasic copolymer)
C10539	Protective cap	● PP (Polypropylene)
C10808	Filter 15 µm drip chamber 21	● Filter: PA 6.6 (Polyamide) = Nylon Filter ● housing: PA 6 (Polyamide),
C10836	ventilation filter 3µm	● Membrane Type: Hydrophobic Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ● Membrane Support Type: Non-woven Polyester
C11504	Solvent / Lubricant	● Isopropyle alcohol
C11511	Lubricant	● Aminofunctional Siloxane
J1010810		
C10405	roller clamp housing, orange	● PS (polystyrene), ● PS (Polystyrene) with silicon slip-masterbatch; ● PS (Polystyrene),
C10412	Roller clamp wheel	● PS (polystyrene); ● 2.5% PS (Polystyrene) with silicone

Article code	Article name	Material
		slipmasterbatch; Alternative to C90348 for FKN only: • PS (Polystyrene),
J1010902		
C11503	Solvent	• THF (Tetrahydrofuran), also 1,4-Epoxy butane,
C11515	Solvent	• MEK (Methylethylketone)
C30631	Connection piece	• MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) • (2%) carrier material ABS (Acrylonitrile Butadiene Styrene) with anorganic and organic colour pigments (blue)
C30632	Pump connector upstream	• MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) • (2%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30650	Antikink	• PP (Polypropylene) • (2%) Carrier material PP with organic and inorganic colourants
C31646	Ambix Nova Clamp	• Polycarbonate • 2% Masterbatch blue (Carrier: Polycarbonate)
C60936	noDOP PVC-tube $\varnothing$ 3,0 x $\varnothing$ 5,0 x 195mm	• PVC (Polyvinylchloride) with plasticizer TEHTM
J1010920		
C10505	Protective cap lock	• PP (Polypropylene),
C10605	connector Luer Lock	• MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene), alternativ: C90106 MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)
Packaging		
	Folding box 400 x 220 x 95 mm	• Corrugated cardboard
	Compound foil 70 μ m	• PA (Polyamide) / adhesive layer / PE (Polyethylene)
	PE-Band 90 x 30 x 0.05 mm	• HDPE (high density polyethylene)
	White paper for IFU printing	• 100% wood pulp
	Roll Label 105 x 148 mm	• Paper (blank white) and Adhesive
	PE bag with bottomseal, transparent	• LDPE (low density polyethylene) • Film thickness: 45 μ m
	Stretch film	• LLDPE (linear low density polyethylene)
	Stretch film foil plate	• LDPE (low density polyethylene)
	PE (Polyethylene) film broadband	• PE (Polyethylene), white
	Paper Ambix nova ambulatory Sets	
Sterilization		
	Sterilization gas, Nanchang	• Ethylene oxide (90%) + Carbon dioxide (10%)

Article code	Article name	Material
	Indicator tape for EO sterilization	Process indicator tape for Ethylene Oxide (EO) • sterilization
	Sterilization gas, Poland	• alternative: > 99.9% Ethylene oxide