

Article code	M46430550Y
Article name acc. DoC	MS 31
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C11503	Solvent	THF (Tetrahydrofuran), 1,4-Epoxy butane,
C30905	Tube 2.94 x 4.09 mm	PVC (Polyvinylchloride), soft, noDOP, with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
C31514	Spring clamp	●Housing: POM (Polyoxymethylene) ●Colourant ●Spring: stainless steel
J1010110		
C10178	Drip chamber upper part 21 E	●polystyrene - Colourant white, 50 % TiO₂ in HIPS (high impact polystyrene)
C10221	moulded ring white for drip chamber 21	●PS (Polystyrene) ●Colourant
C10265	Drip chamber lower part 21 V92	SBS (Styrene Butadiene Copolymer)
C10345	air ventilation flap	●PP (Polypropylene) ●propene copolymer (heterophasic copolymer)
C10539	Protective cap	PP (Polypropylene)
C10808	Filter 15 µm drip chamber 21	●Filter: PA 6.6 (Polyamide) = Nylon Filter ●housing: PA 6 (Polyamide)
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ●Membrane Support Type: Non-woven Polyester
C11504	Solvent / Lubricant	●Isopropyle alcohol
C11511	Lubricant	●Aminofunctional Siloxane
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch; ●PS (Polystyrene)
C10412	Roller clamp wheel	●PS (polystyrene); ●PS (Polystyrene) with silicone slipmasterbatch; ●PS (Polystyrene)
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane

Article code	Article name	Material
C30605	Revolving Luer Lock	●luer cone: PVC (Polyvinylchloride) ●screwcap: ABS (Acrylonitrile Butadiene Styrene)
J1040246		
C31536	Silicone fluid MED-420-350 cP	Liquid copolymer
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
	White paper for IFU printing	●100% wood pulp
	Roll Label 105 x 148 mm	●Paper (blank white) and Adhesive
	Adhesive tape	●BOPP (biaxially oriented polypropylene) film and rubber adhesive
	Carton 400 x 265 x 205 mm	●Cardboard
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Paper for M46430550Y	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide