

Article code	M46441920Y
Article name acc. DoC	VL ST02 AirBlock
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10478	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C30883	Roller clamp safe	●Housing: PS (Polystyrene), colourant white ●Wheel: ABS (Acrylonitrile Butadiene Styrene), colourant blue
C30945	Tube 3.0 x 4.1 mm	●PVC (polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30945S	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
J1010486		
C10221	moulded ring white for drip chamber 21	●PS (Polystyrene) ●Colourant, carrier material PS (Polystyrene), contains 50% TiO₂
C10345	air ventilation flap	●1-propene homopolymer: PP (Polypropylene) ●1-propene copolymer (heterophasic copolymer)
C10539	Protective cap	●PP (Polypropylene)
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ●Membrane Support Type: Non-woven Polyester
C30879	Drip chamber upper part 21 safe	●PS (Polystyrene) ●carrier material S containing anorganic and organic colourants
C30881	Drip chamber lower part 21 safe	●SBS (Styrene Butadiene Copolymer)
C30882	Membrane safe 4 µm	●PES (Polyethersulfone)
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30605	Revolving Luer Lock	●luer cone: PVC (Polyvinylchloride) ●screwcap: ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)

Article code	Article name	Material
C11504	Solvent / Lubricant	●Isopropyle alcohol
C30618	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	●PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	●Raumedic-SIK 8806, Shore A 50-55
J1040246		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●colour batch light blue
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	White paper for IFU printing	●100% wood pulp
	Adhesive tape	●BOPP (biaxially oriented polypropylene) film and rubber adhesive
	Carton 400 x 265 x 205 mm	●Cardboard
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Paper for M46441920Y	●Paper
	PE bag	●LDPE (Low density polyethylene)
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide