

Article code	M46442200Y
Article name acc. DoC	VL PA02 Paediatric
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10478	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C30615	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30945	Tube 3.0 x 4.1 mm	●PVC (polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30945S	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C30949	Tube 2.5 x 4.22 mm	●PVC (Polyvinylchloride)
C41133	Drip chamber complete	●Spike: ABS (Acrylonitrile Butadiene Styrene) ●Protective cap: LDPE (low density polyethylene) ●Drip chamber: PVC (Polyvinylchloride) ●Micro dripper: PC (Polycarbonat) + Stainless steel ●Disk filter: ABS (Acrylonitrile Butadiene Styrene) + PA (Polyami
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene); ●PS (Polystyrene) with silicone slip-masterbatch; ●masterbatch orange; - Alternative to C90348 for FKN only: ●PS (Polystyrene)
C10412	Roller clamp wheel	●PS (polystyrene); ●2.5% PS (Polystyrene) with silicone slip-masterbatch; - Alternative to C90348 for FKN only: ●PS (Polystyrene)
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane

Article code	Article name	Material
C30605	Revolving Luer Lock	•luer cone: PVC (Polyvinylchloride) •screwcap: ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	•Isopropyle alcohol
C30618	Tube adaptor	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) • (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C30674	Tube adaptor	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) • (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	•PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	•Raumedic-SIK 8806, Shore A 50-55
J1040246		
C31536	Silicone fluid MED-420-350 cP	•Liquid copolymer
C31636	K-Zero Female Luer	•Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) •colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	•Liquid Silicone Rubber •Phenyl Oil
Packaging		
	Tubular foil for packaging	•Film 500m x 350 mm x 40 µm •Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Adhesive tape paper	•gummed paper tapes - Polyeste filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	•HDPE (high density polyethylene)
	Label sticker, size: 104 x 148 mm	•Paper for laserprinting
	White paper for IFU printing	•100% wood pulp
	Carton 400 x 265 x 205 mm	•Cardboard
	Compound foil 90µm	•Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)

Article code	Article name	Material
	Paper for M46442200Y	●Paper
Sterilization		
	Sterilization gas, Nanchang	●Ethylene oxide (90%) + Carbon dioxide (10%)
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization