

Article code	M46442500Y
Article name acc. DoC	Infusion Set VL ST 22
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10463	Roberts clamp, white	• ABS (Acrylonitrile Butadiene Styrene)
C10478	Lock plate	• POM (Polyoxymethylene) • carrier material Ionomer with anorganic and organic colour pigments (blue)
C10478B	Lock plate	• POM (Polyoxymethylene) • carrier material Ionomer with anorganic and organic colour pigments (blue)
C11503	Solvent	• THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	• MEK (Methylethylketone)
C30945	Tube 3.0 x 4.1 mm	• PVC (polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30945S	Tube 3.0 x 4.1 mm	• PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C31580	solvent mixture	• 61,3 % THF (Tetrahydrofuran) • 37,2 % MEK (Methylethylketone) • 1,5 % MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
J1010110		
C10178	Drip chamber upper part 21 E	• PS (polystyrene) • Colourant white, 50 % TiO₂ in HIPS (high impact polystyrene)
C10221	moulded ring white for drip chamber 21	• (98%) PS (Polystyrene) • (2%) Colourant, carrier material PS (Polystyrene), contains 50% TiO₂
C10265	Drip chamber lower part 21 V92	• SBS (Styrene Butadiene Copolymer)
C10345	air ventilation flap	• (50%)propene homopolymer: PP (Polypropylene) • (50%)propene copolymer (heterophasic copolymer)
C10539	Protective cap	• PP (Polypropylene)
C10808	Filter 15 µm drip chamber 21	• Filter: PA 6.6 (Polyamide) = Nylon Filter • housing: PA 6 (Polyamide)
C10836	ventilation filter 3µm	• Membrane Type: Hydrophobic • Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane • Membrane Support Type: Non-woven Polyester
C11504	Solvent / Lubricant	• Isopropyle alcohol
C11511	Lubricant	• Aminofunctional Siloxane
J1010810		
C10405	roller clamp housing, orange	• PS (polystyrene)PS (Polystyrene) with silicone slip-

Article code	Article name	Material
		masterbatch • masterbatch orange
C10412	Roller clamp wheel	• PS (polystyrene) • 2.5% PS (Polystyrene) with silicone slip-masterbatch
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	• Housing: PP (Polypropylene) + Masterbatch • Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30605	Revolving Luer Lock	• luer cone: PVC (Polyvinylchloride) • screwcap: ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	• MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	• Isopropyle alcohol
C30618	Tube adaptor	• MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) • (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	• MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	• MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) • (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes(light green)
C30674	Safety ring	• PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	• Raumedic-SIK 8806, Shore A 50-55
J1040246		
C31536	Silicone fluid MED-420-350 cP	• Liquid copolymer
C31636	K-Zero Female Luer	• Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	• MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) • colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	• Liquid Silicone Rubber • Phenyl Oil
J1040330		•
C10564	Protective cap LL pos. with 1.2 µm filter	• Housing: PP (Polypropylene) +Masterbatch • Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C11503	Solvent	• THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	• MEK (Methylethylketone)
		• luer cone: PVC (Polyvinylchloride)

Article code	Article name	Material
C30605	Revolving Luer Lock	● screwcap: ABS (Acrylonitrile ButadieneStyrene)
C30945	Tube 3.0 x 4.1 mm	● PVC (polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30945S	Tube 3.0 x 4.1 mm	● PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C31536	Silicone fluid MED-420-350 cP	● Liquid copolymer
C31592	Tinopal Optical Brightener	● 2,5-thiophenediylbis(5-tert-butyl-1,3-benzoxazole), color: Slight greenish yellow
C31636	K-Zero Female Luer	● Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	● Liquid Silicone Rubber ● Phenyl Oil
Packaging		
	Tubular foil for packaging	● Film 500m x 350 mm x 40 µm ● Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	PE-Band 90 x 30 x 0.05 mm	● HDPE (high density polyethylene)
	Label sticker, size: 104 x 148 mm	● Paper for laserprinting
	White paper for printing, IFU	● 100% wood pulp
	Adhesive tape	● BOPP (biaxially oriented polypropylene) film and rubber adhesive
	Carton 400 x 265 x 205 mm	● Cardboard
	Stretch film	● LLDPE (linear low-density polyethylene)
	Stretch film foil plate	● LDPE (low density polyethylene)
	Compound foil 90µm	● Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Paper for M46441980Y, 133,67x267 mm	● Paper
Sterilization		
	Sterilization gas	● Ethylene oxide
	Indicator tape for EO sterilization	● Process indicator tape for Ethylene Oxide (EO) sterilization