

Article code	M46442800S
Article name acc. DoC	VL TR00 Transfusion
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10478	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C10478B	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C30945S	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C31580	solvent mixture	●61,3 % THF (Tetrahydrofuran) ●37,2 % MEK (Methylethylketone) ●1,5 % MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
J1010446		
C10104	Drip chamber upper part blood bag	●PS (Polystyrene)
	Drip chamber lower part 18	●SBS (Styrene Butadiene Copolymer)
C10216	Moulded ring white for drip chamber 18	●(98%) PS (Polystyrene) ●2% Colourant, carrier material PS (Polystyrene), contains 50% TiO2
C10539	Protective cap	●PP (Polypropylene)
C10810	Tube filter 200 µm	●Filter housing: Resin PS (Polystyrene) ●Filter: Mesh PA 6.6 (Polyamide) = Nylon
C11504	Solvent / Lubricant	●Isopropyle alcohol
	Lubricant	●Aminofunctional Siloxane
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch; ●masterbatch orange; - Alternative tofor FKN only: ●PS (Polystyrene)

Article code	Article name	Material
C10412	Roller clamp wheel	●PS (polystyrene); ●2.5% PS (Polystyrene) with silicone slip-masterbatch; - Alternative to for FKN only: ●PS (Polystyrene)
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
	Revolving Luer Lock	●luer cone: PVC (Polyvinylchloride) ●screwcap: ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	●T MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C30618	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●(3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●(3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	●PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	●Raumedic-SIK 8806, Shore A 50-55
Packaging		
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
C11521	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
C11600	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
C11945	White paper for IFU printing	●100% wood pulp
C12184	Adhesive tape	●BOPP (biaxially oriented polypropylene) film and rubber adhesive

Article code	Article name	Material
C31016	PE bag with bottomseal, transparent	●LDPE (low density polyethylene) ●Film thickness: 45 µm
C31025A	Carton 400 x 265 x 205 mm	●Cardboard
	Stretch film	●LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	●LDPE (low density polyethylene)
C31055	PE (Polyethylene) film broadband	●PE (Polyethylene), white
C31301	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
C62389CN	Paper "Volumat Line Transfusion"	●Paper
C62389PL	Paper Volumat Line	●Paper
Sterilization		
C11520CN	Sterilization gas, Nanchang	●Ethylene oxide (90%) + Carbon dioxide (10%)
●THF (Tetrahydrofuran), also 1,4-Epoxy butane	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
C31500	Sterilization gas, Poland	●100% Ethylene oxide, liquefied grade one (GAZ 4-134) ●alternative: > 99.9% Ethylene oxide