

Article code	M46442800S
Article name acc. DoC	Transfusion Set VL TR 00
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10478	Lock plate	•POM (Polyoxymethylene) •carrier material lonomer with anorganic and organic colour pigments (blue)
C10478B	Lock plate	•POM (Polyoxymethylene) •carrier material lonomer with anorganic and organic colour pigments (blue)
C11503	Solvent	•THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	•MEK (Methylethylketone)
C30945S	Tube 3.0 x 4.1 mm	•PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C31580	solvent mixture	•61,3 % THF (Tetrahydrofuran) •37,2 % MEK (Methylethylketone) •1,5 % MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
J1010446		
C10104	Drip chamber upper part blood bag	•PS (Polystyrene)
C10203	Drip chamber lower part 18	•SBS (Styrene Butadiene Copolymer)
C10216	Moulded ring white for drip chamber 18	• (98%) PS (Polystyrene) •2% Colourant, carrier material PS (Polystyrene), contains 50% TiO₂
C10539	Protective cap	•PP (Polypropylene)
C10810	Tube filter 200 µm	•Filter housing: Resin PS (Polystyrene) •Filter: Mesh PA 6.6 (Polyamide) = Nylon
C11504	Solvent / Lubricant	•Isopropyle alcohol
C11511	Lubricant	•Aminofunctional Siloxane
J1010810		
C10405	roller clamp housing, orange	•PS (polystyrene) •PS (Polystyrene) with silicone slip- masterbatch •masterbatch orange
C10412	Roller clamp wheel	•PS (polystyrene) •2.5% PS (Polystyrene) with silicone slip-masterbatch
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	•Housing: PP (Polypropylene) + Masterbatch •Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane

Article code	Article name	Material
C30605	Revolving Luer Lock	•luer cone: PVC (Polyvinylchloride) •screwcap: ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	•Isopropyle alcohol
C30618	Tube adaptor	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) • (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	•PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	•Raumedic-SIK 8806, Shore A 50-55
Packaging		
	Tubular foil for packaging	•Film 500m x 350 mm x 40 µm •Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	PE-Band 90 x 30 x 0.05 mm	•HDPE (high density polyethylene)
	Label sticker, size: 104 x 148 mm	•Paper for laserprinting
	White paper for printing, IFU	100% wood pulp
	Adhesive tape	•BOPP (biaxially oriented polypropylene) film and rubber adhesive
	PE bag with bottomseal, transparent	•LDPE (low density polyethylene) •Film thickness: 45 µm
	Carton 400 x 265 x 205 mm	•Cardboard
	Stretch film	•LLDPE (linear low-density polyethylene)
	Stretch film foil plate	•LDPE (low density polyethylene)
	PE (Polyethylene) film broadband	•PE (Polyethylene), white
	Compound foil 90µm	•Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Paper "Volumat Line Transfusion"	•Paper
	Paper Volumat Line	•Paper
Sterilization		
	Sterilization gas	•Ethylene oxide
	Indicator tape for EO sterilization	•Process indicator tape for Ethylene Oxide (EO) sterilization