

Article code	M46442900Y
Article name acc. DoC	Transfusion Set VL TR 00
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10463	Roberts clamp, white	ABS (Acrylonitrile Butadiene Styrene)
C10478	Lock plate	• POM (Polyoxymethylene) • carrier material Ionomer with anorganic and organic colour pigments (blue)
C11503	Solvent	• THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	MEK (Methylethylketone)
C30706	Back check valve	• Solvent: TEF (Tetrahydrofuran) • Check Valve: MABS (Methymethacrylate Acrylonitrile Butadiene Styrene), 2% Masterbatch Valve membrane: Silicone
C30945	Tube 3.0 x 4.1 mm	• PVC (polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30945S	Tube 3.0 x 4.1 mm	• PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C31580	solvent mixture	• 61,3 % THF (Tetrahydrofuran) • 37,2 % MEK (Methylethylketone) • 1,5 % MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
J1010299		
C10178	Drip chamber upper part 21 E	• PS (polystyrene) • Colourant white, 50 % TiO₂ in HIPS (high impact polystyrene)
C10221	moulded ring white for drip chamber 21	• (98%) PS (Polystyrene) • (2%) Colourant, carrier material PS (Polystyrene), contains 50% TiO₂
C10243	Lower part 21 f.blood filter	SBS (Styrene-Butadiene-Copolymer)
C10345	air ventilation flap	• (50%)propene homopolymer: PP (Polypropylene) • (50%)propene copolymer (heterophasic copolymer)
C10539	Protective cap	PP (Polypropylene)
C10813	Tube filter 200 µm	• Filter housing: ABS (Acrylonitrile Butadiene Styrene) Filter cloth: PA 6.6 (Polyamide) = Nylon
C10836	ventilation filter 3µm	• Membrane Type: Hydrophobic • Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane

Article code	Article name	Material
		Membrane Support Type: Non-woven Polyester
C11504	Solvent / Lubricant	Isopropyle alcohol
C11511	Lubricant	Aminofunctional Siloxane
J1010810		
C10405	roller clamp housing, orange	•PS (polystyrene) •PS (Polystyrene) with silicone slip- masterbatch •masterbatch orange
C10412	Roller clamp wheel	•PS (polystyrene) •2.5% PS (Polystyrene) with silicone slip-masterbatch
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	•Housing: PP (Polypropylene) + Masterbatch •Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30605	Revolving Luer Lock	•luer cone: PVC (Polyvinylchloride) •screwcap: ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	•MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	•Isopropyle alcohol
C30618	Tube adaptor	•MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) • (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	•MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	•MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	•PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	•Raumedic-SIK 8806, Shore A 50-55
J1040246		
C31536	Silicone fluid MED-420-350 cP	Liquid copolymer
C31636	K-Zero Female Luer	Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	•MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) - colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	•Liquid Silicone Rubber - Phenyl Oil
J1040330		

Article code	Article name	Material
C10564	Protective cap LL pos. with 1.2 µm filter	- Housing: PP (Polypropylene) + Masterbatch - Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C11503	Solvent	THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	MEK (Methylethylketone)
C30605	Revolving Luer Lock	- luer cone: PVC (Polyvinylchloride) - screwcap: ABS (Acrylonitrile Butadiene Styrene)
C30945	Tube 3.0 x 4.1 mm	PVC (polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30945S	Tube 3.0 x 4.1 mm	PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C31536	Silicone fluid MED-420-350 cP	Liquid copolymer
C31592	Tinopal Optical Brightener	2,5-thiophenediylbis(5-tert-butyl-1,3-benzoxazole), color: Slight greenish yellow
C31636	K-Zero Female Luer	Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	- MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) - colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	- Liquid Silicone Rubber - Phenyl Oil
Packaging		
	Tubular foil for packaging	Film 500m x 350 mm x 40 µm Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Adhesive tape paper	- gummed paper tapes - Polyeste filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	HDPE (high density polyethylene)
	White paper for printing, IFU	100% wood pulp
	Roll Label 105 x 148 mm	Paper (blank white) and Adhesive
	Carton 400 x 265 x 205 mm	Cardboard
	Compound foil 90µm	Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Paper for M46442900Y	Paper
Sterilization		
	Sterilization gas	Ethylene oxide
	Indicator tape for EO sterilization	Process indicator tape for Ethylene Oxide (EO) sterilization