

Material List - M46443100S

Article code	M46443100S
Article name acc. DoC	Transfusion Set VL SP 22
Class acc. Directive 93/42/EEC	Ila

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10478	Lock plate	●POM (Polyoxymethylene) ● carrier material Ionomer with anorganic and organic colour pigments (blue)
C10478B	Lock plate	●POM (Polyoxymethylene) ● carrier material Ionomer with anorganic and organic colour pigments (blue)
C10702	3-way piece	● MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C30727	needle free access device	●Housing: PC (Polycarbonate), Colourant and transparent part ●Stopper: Silicone
C30945	Tube 3.0 x 4.1 mm	●PVC (polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C31580	solvent mixture	●61,3 % THF (Tetrahydrofuran) ●37,2 % MEK (Methylethylketone) ●1,5 % MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
J1010465US		
C10103	Drip chamber upper part	●PS (Polystyrene)
C10203	Drip chamber lower part 18	●SBS (Styrene Butadiene Copolymer)
	Moulded ring white for drip chamber 18	●(98%) PS (Polystyrene) ●2% Colourant, carrier material PS (Polystyrene), contains 50% TiO ₂
C10837	Tube filter 200 µm	●Housing: Polyamide 6 Akulon K222-D ●Filter: Mesh PA 6.6 (Polyamide 6.6)
C11504	Solvent / Lubricant	●Isopropyle alcohol
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene) ●PS (Polystyrene) with silicone

Material List - M46443100S

Article code	Article name	Material
		●masterbatch orange; Alternative to C90348 for FKN only: ●PS (Polystyrene)
C10412	Roller clamp wheel	●PS (polystyrene); ●2.5% PS (Polystyrene) with silicone slip-masterbatch; Alternative to for FKN only: ●PS (Polystyrene)
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30605	Revolving Luer Lock	●luer cone: PVC (Polyvinylchloride) ●screwcap: ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene),
C11504	Solvent / Lubricant	●Isopropyle alcohol
	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●(3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●(3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	●PP (Polypropylene)
	SIK-Tube Ø3,41xØ5,11	●Raumedic-SIK 8806, Shore A 50-55
J1020091		
C10539	Protective cap	●PP (Polypropylene)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C11511	Lubricant	●Aminofunctional Siloxane
C21791	Vented spike body	●PS (polystyrene) ●Colourant white, 50 % TiO₂ in HIPS (high impact polystyrene)
Packaging		

Material List - M46443100S

Article code	Article name	Material
	Tubular foil for packaging	● Film 500m x 350 mm x 40 µm ● Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
C11517	Adhesive tape paper	● gummed paper tapes ● Polyester filament yarn for strengthening
C11556	PE banderole 110 x 180 x 0.05 mm	● LDPE (low density polyethylene)
C11945	White paper for IFU printing	● 100% wood pulp
C12182	Roll Label 105 x 148 mm	● Paper (blank white) and Adhesive
C31016	PE bag with bottomseal, transparent	● LDPE (low density polyethylene) ● Film thickness: 45 µm
C31025A	Carton 400 x 265 x 205 mm	● Cardboard
C31045	Stretch film	● LLDPE (linear low density polyethylene)
● THF (Tetrahydrofuran), also 1,4-Epoxy butane	Stretch film foil plate	● LDPE (low density polyethylene)
C31055	PE (Polyethylene) film broadband	● PE (Polyethylene), white
C31301	Compound foil 90µm	● Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
C62389	Paper Transfusion Set 133 x 267 mm	● Paper
C64576	Paper Volumat Line	● paper
Sterilization		
C11520CN	Sterilization gas, Nanchang	● Ethylene oxide (90%) + Carbon dioxide (10%)
C11525	Indicator tape for EO sterilization	● Process indicator tape for Ethylene Oxide (EO) sterilization
C31500	Sterilization gas, Poland	● 100% Ethylene oxide, liquefied grade one (GAZ 4-134) ● alternative: > 99.9% Ethylene oxide