

Article code	M46444061Y
Article name acc. DoC	VL ON42 SC PUR material, 3-way-stopcock
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10478US	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C10511	protective cap LL	●PS (polystyrene) with butadien HIPS (high impact polystyrene) ●PS (polystyrene) with silicone
C10605	connector Luer Lock	●MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)
C10942	Tube 3.0 x 4.1 mm	●Polyether Polyurethane
C10942A	Polyurethane Tube Ø3,00 x Ø4,1 mm	●Polyether Polyurethane
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C11519	Solvent	●DCM (Dichlormethane) + 0.2 % Ethanol
C11553	Solvent	●CHN (Cyclohexanone)
C30756	Three-Way-Stopcock	● 3-arm piece (core) PE (Polyethylene, high density) + colourant ● 4-way body PC (Polycarbonate) ● Colourant blue ● Lubricant Silicone oil
C30911	Tube 2.94 x 4.02 mm	●PP (Polypropylene)
J1010422		
C10148	Drip chamber upper part	●PS (Polystyrene) ●Colourant white, 50 % TiO ₂ in HIPS (high impact polystyrene) optional ●PS (polystyrene) with silicone
C10213	Drip chamber lower part 21 long	●SBS (Styrene-Butadiene-Copolymer)
C10221	moulded ring white for drip chamber 21	●PS (Polystyrene) ●Colourant, carrier material PS (Polystyrene), contains 50% TiO ₂
C10808	Filter 15 µm drip chamber 21	●Filter: PA 6.6 (Polyamide) = Nylon Filter ●housing: PA 6 (Polyamide)
J1010810		

Article code	Article name	Material
C10405	roller clamp housing, orange	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch; ●masterbatch orange
C10412	Roller clamp wheel	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch
J1010915		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30675	Revolving Luer Lock	●ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C30618	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	●PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	●Raumedic-SIK 8806, Shore A 50-55
J1020091		
C10345	air ventilation flap	●1-propene homopolymer: PP (Polypropylene) ●1-propene copolymer (heterophasic copolymer)
C10539	Protective cap	●PP (Polypropylene)
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ●Membrane Support Type: Non-woven Polyester
C11504	Solvent / Lubricant	●Isopropyle alcohol
C11511	Lubricant	●Aminofunctional Siloxane
C21791	Vented spike body	●PS (polystyrene) ●Colourant white, 50 % TiO₂ in HIPS (high impact polystyrene)
J1040246		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer

Article code	Article name	Material
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●colour batch light blue
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Adhesive tape paper	●gummed paper tapes Polyeste filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
	White paper for IFU printing	●100% wood pulp
	Roll Label 105 x 148 mm	●Paper (blank white) and Adhesive
	Carton 400 x 265 x 205 mm	●Cardboard
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Paper for M46444061Y	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide