

|                                |                                                             |
|--------------------------------|-------------------------------------------------------------|
| <b>Article code</b>            | M46444071Y                                                  |
| <b>Article name acc. DoC</b>   | VL ON42 OP SC Light-protected, PUR material, 3-way-stopcock |
| <b>Class acc. MDR 2017/745</b> | Ila                                                         |

| <b>Article code</b> | <b>Article name</b>     | <b>Material</b>                                                                                                                                                                       |
|---------------------|-------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| C10463              | Roberts clamp, white    | ●ABS (Acrylonitrile Butadiene Styrene)                                                                                                                                                |
| C10478US            | Lock plate              | ●POM (Polyoxymethylene)<br>●carrier material Ionomer with anorganic and organic colour pigments (blue)                                                                                |
| C10511              | protective cap LL       | ●PS (polystyrene) with butadien HIPS (high impact polystyrene)<br>●PS (polystyrene) with silicone                                                                                     |
| C11503              | Solvent                 | ●THF (Tetrahydrofuran), also 1,4-Epoxy butane                                                                                                                                         |
| C11515              | Solvent                 | ●MEK (Methylethylketone)                                                                                                                                                              |
| C11519              | Solvent                 | ●DCM (Dichlormethane) + 0.2 % Ethanol                                                                                                                                                 |
| C11553              | Solvent                 | ●CHN (Cyclohexanone)                                                                                                                                                                  |
| C30629              | Tube adaptor            | ●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)<br>●carrier material MABS (Terlux) with organic colour pigments (yellow)                                                    |
| C30633              | Male connector LL       | ●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)<br>●carrier material MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)(Terlux) with organic colour pigments (yellow) |
| C30756              | Three-Way-Stopcock      | ● 3-arm piece (core) PE (Polyethylene, high density) + colourant<br>● 4-way body PC (Polycarbonate)<br>● Colourant blue<br>● Lubricant Silicone oil                                   |
| C30911              | Tube 2.94 x 4.02 mm     | ●PP (Polypropylene)                                                                                                                                                                   |
| C30947              | Tube 3 x 4.1 mm         | ●PUR (Polyurethane) + 3,85% yellow UV stabilized Masterbatch PUR yellow                                                                                                               |
| J1010484            |                         |                                                                                                                                                                                       |
| C10148              | Drip chamber upper part | ●PS (Polystyrene)<br>●Colourant white, 50 % TiO <sub>2</sub> in HIPS (high impact polystyrene)<br>optional<br>●PS (polystyrene) with silicone                                         |

| Article code | Article name                              | Material                                                                                                                                                                                                            |
|--------------|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| C10221       | moulded ring white for drip chamber 21    | <ul style="list-style-type: none"> <li>●PS (Polystyrene)</li> <li>●Colourant, carrier material PS (Polystyrene), contains 50% TiO<sub>2</sub></li> </ul>                                                            |
| C10274       | Drip chamber lower part                   | <ul style="list-style-type: none"> <li>●SBS (Styrene-Butadiene-Copolymer)</li> <li>●masterbatch with UV-stabilizer Carrier: PS (Polystyrene)</li> </ul>                                                             |
| C10808       | Filter 15 µm drip chamber 21              | <ul style="list-style-type: none"> <li>●Filter: PA 6.6 (Polyamide) = Nylon Filter</li> <li>●housing: PA 6 (Polyamide)</li> </ul>                                                                                    |
| J1010810     |                                           |                                                                                                                                                                                                                     |
| C10405       | roller clamp housing, orange              | <ul style="list-style-type: none"> <li>●PS (polystyrene)</li> <li>●PS (Polystyrene) with silicone slip-masterbatch;</li> <li>●masterbatch orange</li> </ul>                                                         |
| C10412       | Roller clamp wheel                        | <ul style="list-style-type: none"> <li>●PS (polystyrene);</li> <li>●PS (Polystyrene) with silicone slip-masterbatch</li> </ul>                                                                                      |
| J1010915     |                                           |                                                                                                                                                                                                                     |
| C10564       | Protective cap LL pos. with 1.2 µm filter | <ul style="list-style-type: none"> <li>●Housing: PP (Polypropylene) + Masterbatch</li> <li>●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane</li> </ul>            |
| C30675       | Revolving Luer Lock                       | <ul style="list-style-type: none"> <li>●ABS (Acrylonitrile Butadiene Styrene)</li> </ul>                                                                                                                            |
| J1010938     |                                           |                                                                                                                                                                                                                     |
| C10477       | Rotation clamp base                       | <ul style="list-style-type: none"> <li>●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)</li> </ul>                                                                                                         |
| C11504       | Solvent / Lubricant                       | <ul style="list-style-type: none"> <li>●Isopropyle alcohol</li> </ul>                                                                                                                                               |
| C30618       | Tube adaptor                              | <ul style="list-style-type: none"> <li>●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)</li> <li>●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)</li> </ul> |
| C30638       | Tube adaptor                              | <ul style="list-style-type: none"> <li>●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)</li> </ul>                                                                                                        |
| C30673       | Tube adaptor                              | <ul style="list-style-type: none"> <li>●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)</li> <li>●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)</li> </ul> |
| C30674       | Safety ring                               | <ul style="list-style-type: none"> <li>●PP (Polypropylene)</li> </ul>                                                                                                                                               |
| C30946       | SIK-Tube Ø3,41xØ5,11                      | <ul style="list-style-type: none"> <li>●Raumedic-SIK 8806, Shore A 50-55</li> </ul>                                                                                                                                 |
| J1020091     |                                           |                                                                                                                                                                                                                     |
| C10345       | air ventilation flap                      | <ul style="list-style-type: none"> <li>●1-propene homopolymer: PP (Polypropylene)</li> <li>●1-propene copolymer (heterophasic copolymer)</li> </ul>                                                                 |
| C10539       | Protective cap                            | <ul style="list-style-type: none"> <li>●PP (Polypropylene)</li> </ul>                                                                                                                                               |

| Article code  | Article name                        | Material                                                                                                                                                                                                     |
|---------------|-------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| C10836        | ventilation filter 3µm              | <ul style="list-style-type: none"> <li>●Membrane Type: Hydrophobic</li> <li>●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane</li> <li>●Membrane Support Type: Non-woven Polyester</li> </ul> |
| C11504        | Solvent / Lubricant                 | ●Isopropyle alcohol                                                                                                                                                                                          |
| C11511        | Lubricant                           | ●Aminofunctional Siloxane                                                                                                                                                                                    |
| C21791        | Vented spike body                   | <ul style="list-style-type: none"> <li>●PS (polystyrene)</li> <li>●Colourant white, 50 % TiO<sub>2</sub> i n HIPS (high impact polystyrene)</li> </ul>                                                       |
| J1040246      |                                     |                                                                                                                                                                                                              |
| C31536        | Silicone fluid MED-420-350 cP       | ●Liquid copolymer                                                                                                                                                                                            |
| C31636        | K-Zero Female Luer                  | ●Eastman Tritan Copolyester MX731                                                                                                                                                                            |
| C31638        | K-Zero Y-Site Body                  | <ul style="list-style-type: none"> <li>●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)</li> <li>●colour batch light blue</li> </ul>                                                               |
| C31639        | K-Zero Microbial Barrier            | <ul style="list-style-type: none"> <li>●Liquid Silicone Rubber</li> <li>●Phenyl Oil</li> </ul>                                                                                                               |
| Packaging     |                                     |                                                                                                                                                                                                              |
|               | Tubular foil for packaging          | <ul style="list-style-type: none"> <li>●Film 500m x 350 mm x 40 µm</li> <li>●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent</li> </ul>                   |
|               | Adhesive tape paper                 | <ul style="list-style-type: none"> <li>●gummed paper tapes</li> <li>Polyeste filament yarn for strengthening</li> </ul>                                                                                      |
|               | PE-Band 90 x 30 x 0.05 mm           | ●HDPE (high density polyethylene)                                                                                                                                                                            |
|               | White paper for IFU printing        | ●100% wood pulp                                                                                                                                                                                              |
|               | Roll Label 105 x 148 mm             | ●Paper (blank white) and Adhesive                                                                                                                                                                            |
|               | Carton 400 x 265 x 205 mm           | ●Cardboard                                                                                                                                                                                                   |
|               | Stretch film                        | ●LLDPE (linear low density polyethylene)                                                                                                                                                                     |
|               | Stretch film foil plate             | ●LDPE (low density polyethylene)                                                                                                                                                                             |
|               | Compound foil 90µm                  | ●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)                                                                                                                                     |
|               | Paper for M46444071Y                | ●Paper                                                                                                                                                                                                       |
| Sterilization |                                     |                                                                                                                                                                                                              |
|               | Indicator tape for EO sterilization | ●Process indicator tape for Ethylene Oxide (EO) sterilization                                                                                                                                                |
|               | Sterilization gas                   | ●Ethylene oxide                                                                                                                                                                                              |