

|                                       |                             |
|---------------------------------------|-----------------------------|
| <b>Article code</b>                   | M46444080S                  |
| <b>Article name acc. DoC</b>          | Infusion Set VL ON 42 OP NF |
| <b>Class acc. Directive 93/42/EEC</b> | IIa                         |

| Article code | Article name                              | Material                                                                                                                                         |
|--------------|-------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| C10463       | Roberts clamp, white                      | ●ABS (Acrylonitrile Butadiene Styrene)                                                                                                           |
| C10478US     | Lock plate                                | ●POM (Polyoxymethylene)<br>●carrier material Ionomer with anorganic and organic colour pigments (blue)                                           |
| C11503       | Solvent                                   | ●THF (Tetrahydrofuran), also 1,4-Epoxy butane                                                                                                    |
| C11515       | Solvent                                   | ●MEK (Methylethylketone)                                                                                                                         |
| C11553       | Solvent                                   | ●CHN (Cyclohexanone)                                                                                                                             |
| C30629       | Tube adaptor                              | ●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)<br>●carrier material MABS (Terlux) with organic colour pigments (yellow)               |
| C30727       | needle free access device                 | ●Housing: PC (Polycarbonate), Colourant and transparent part<br>●Stopper: Silicone                                                               |
| C30911       | Tube 2.94 x 4.02 mm                       | ●PP (Polypropylene)                                                                                                                              |
| C30947       | Tube 3 x 4.1 mm                           | ●PUR (Polyurethane) + yellow UV stabilized Masterbatch PUR yellow                                                                                |
| J1010204     |                                           |                                                                                                                                                  |
| C10148       | Drip chamber upper part                   | ●PS (Polystyrene)<br>●Colourant white, 50 % TiO <sub>2</sub> in HIPS (high impact polystyrene)<br>optional<br>●PS (polystyrene) with silicone    |
| C10221       | moulded ring white for drip chamber 21    | ●PS (Polystyrene)<br>●Colourant, carrier material PS (Polystyrene), contains 50% TiO <sub>2</sub>                                                |
| C10266       | Drip chamber lower part 21 V92, yellow    | ●SBS (Styrene Butadiene Copolymer)<br>●masterbatch with UV-stabilizer, yellow, Carrier: PS (Polystyrene)                                         |
| C10808       | Filter 15 µm drip chamber 21              | ●Filter: PA 6.6 (Polyamide) = Nylon Filter<br>●housing: PA 6 (Polyamide)                                                                         |
| J1010810     |                                           |                                                                                                                                                  |
| C10405       | roller clamp housing, orange              | ●PS (polystyrene)<br>●PS (Polystyrene) with silicone slip-masterbatch;<br>●masterbatch orange;                                                   |
| C10412       | Roller clamp wheel                        | ●PS (polystyrene);<br>●PS (Polystyrene) with silicone slip-masterbatch;                                                                          |
| J1010915     |                                           |                                                                                                                                                  |
| C10564       | Protective cap LL pos. with 1.2 µm filter | ●Housing: PP (Polypropylene) + Masterbatch<br>●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane |
| C30675       | Revolving Luer Lock                       | ●ABS (Acrylonitrile Butadiene Styrene)                                                                                                           |

| Article code  | Article name                 | Material                                                                                                                                                    |
|---------------|------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| J1010938      |                              |                                                                                                                                                             |
| C10477        | Rotation clamp base          | ●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)                                                                                                   |
| C11504        | Solvent / Lubricant          | ●Isopropyle alcohol                                                                                                                                         |
| C30618        | Tube adaptor                 | ●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)<br>●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green) |
| C30638        | Tube adaptor                 | ●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)                                                                                                  |
| C30673        | Tube adaptor                 | ●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)<br>●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green) |
| C30674        | Safety ring                  | ●PP (Polypropylene)                                                                                                                                         |
| C30946        | SIK-Tube Ø3,41xØ5,11         | ●Raumedic-SIK 8806, Shore A 50-55                                                                                                                           |
| J1020091      |                              |                                                                                                                                                             |
| C10345        | air ventilation flap         | ●1-propene homopolymer: PP (Polypropylene)<br>●1-propene copolymer (heterophasic copolymer)                                                                 |
| C10539        | Protective cap               | ●PP (Polypropylene)                                                                                                                                         |
| C10836        | ventilation filter 3µm       | ●Membrane Type: Hydrophobic<br>●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane<br>●Membrane Support Type: Non-woven Polyester              |
| C11504        | Solvent / Lubricant          | ●Isopropyle alcohol                                                                                                                                         |
| C11511        | Lubricant                    | ●Aminofunctional Siloxane                                                                                                                                   |
| C21791        | Vented spike body            | ●PS (polystyrene)<br>●Colourant white, 50 % TiO <sub>2</sub> in HIPS (high impact polystyrene)                                                              |
| Packaging     |                              |                                                                                                                                                             |
|               | PE bag                       | ●LDPE (Low density polyethylene)                                                                                                                            |
|               | Tubular foil for packaging   | ●Film 500m x 350 mm x 40 µm<br>●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent                          |
|               | Adhesive tape paper          | ●gummed paper tapes<br>Polyeste filament yarn for strengthening                                                                                             |
|               | PE-Band 90 x 30 x 0.05 mm    | ●HDPE (high density polyethylene)                                                                                                                           |
|               | White paper for IFU printing | ●100% wood pulp                                                                                                                                             |
|               | Roll Label 105 x 148 mm      | ●Paper (blank white) and Adhesive                                                                                                                           |
|               | Carton 400 x 265 x 205 mm    | ●Cardboard                                                                                                                                                  |
|               | Stretch film                 | ●LLDPE (linear low density polyethylene)                                                                                                                    |
|               | Stretch film foil plate      | ●LDPE (low density polyethylene)                                                                                                                            |
|               | Compound foil 90µm           | ●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)                                                                                    |
|               | Paper Volumat Line           | ●Paper                                                                                                                                                      |
| Sterilization |                              |                                                                                                                                                             |

| Article code | Article name                        | Material                                                       |
|--------------|-------------------------------------|----------------------------------------------------------------|
|              | Sterilization gas                   | ● Ethylene oxide                                               |
|              | Indicator tape for EO sterilization | ● Process indicator tape for Ethylene Oxide (EO) sterilization |