

Article code	M46444230Y
Article name acc. DoC	VL ON22 OP NF Light-protected, PUR material
Class acc. MDR 2017/745	IIa

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10478US	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C30629	Tube adaptor	● MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) ●carrier material MABS (Terlux) with organic colour pigments (yellow)
C30911	Tube 2.94 x 4.02 mm	●PP (Polypropylene)
C30947	Tube 3 x 4.1 mm	●PUR (Polyurethane) + 3,85% yellow UV stabilized Masterbatch PUR yellow
J1010204		
C10148	Drip chamber upper part	● PS (Polystyrene) ● Colourant white, 50 % TiO ₂ in HIPS (high impact polystyrene)optional ●PS (polystyrene) with silicone
C10221	moulded ring white for drip chamber 21	●(98%) PS (Polystyrene) ● (2%) Colourant, carrier material PS (Polystyrene), contains 50% TiO ₂
C10266	Drip chamber lower part 21 V92, yellow	●SBS (Styrene Butadiene Copolymer)(95%) ●masterbatch with UV-stabilizer, 5% yellow, Carrier: PS (Polystyrene) (5%)
C10808	Filter 15 µm drip chamber 21	●Filter: PA 6.6 (Polyamide) = Nylon Filter ●housing: PA 6 (Polyamide)
C11504	Solvent / Lubricant	●Isopropyle alcohol
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch; ● masterbatch orange; ● PS (Polystyrene)
C10412	Roller clamp wheel	●PS (polystyrene); ●2.5% PS (Polystyrene) with silicone slip-masterbatch; ●C90699 PS (Polystyrene)
J1010915		

C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30675	Revolving Luer Lock	●ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C30618	Tube adaptor	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●(3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●(3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	● PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	●Raumedic-SIK 8806, Shore A 50-55
J1020091		
C10539	Protective cap	● PP (Polypropylene)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C11511	Lubricant	●Aminofunctional Siloxane
C21791	Vented spike body	● PS (polystyrene) ●Colourant white, 50 % TiO₂ in HIPS (high impact polystyrene)
J1040246		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31636	K-Zero Female Luer	● Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Adhesive tape paper	●gummed paper tapes Polyeste filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
	White paper for IFU printing	●100% wood pulp

	Roll Label 105 x 148 mm	●Paper (blank white) and Adhesive
	Carton 400 x 265 x 205 mm	●Cardboard
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Paper for M46444230Y, 133,67x267 mm	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	●100% Ethylene oxide, liquefied grade one (GAZ 4-134) ●alternative: > 99.9% Ethylene oxide