

Article code	M46444234
Article name acc. DoC	VL ON22 OP NF AirBlock PUR material, lightprotected
Class acc. MDR 2017/745	IIa

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10478US	Lock plate	● POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C11553	Solvent	●CHN (Cyclohexanone)
C30629	Tube adaptor	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) ● carrier material MABS (Terlux) with organic colour pigments (yellow)
C30883	Roller clamp safe	●Housing: PS (Polystyrene), colourant white ●Wheel: ABS (Acrylonitrile Butadiene Styrene), colourant blue
C30911	Tube 2.94 x 4.02 mm	●PP (Polypropylene)
C30947	Tube 3 x 4.1 mm	●PUR (Polyurethane) + 3,85% yellow UV stabilized Masterbatch PUR yellow
J1010341		
C10193	Tube adaptor, blue	●PS (Polystyrene) ●(3%) carrier material S P containing anorganic and organic colourants
C10221	moulded ring white for drip chamber 21	●(98%) PS (Polystyrene) ● (2%) Colourant, carrier material PS (Polystyrene), contains 50% TiO ₂
C10278	Drip chamber lower part (UT 21 safe)	● SBS (Styrene Butadiene Copolymer) ● masterbatch with UV-stabilizer, 5% yellow, Carrier: PS (Polystyrene) (5%)
C30882	Membrane safe 4 µm	●PES (Polyethersulfone)
J1010915		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30675	Revolving Luer Lock	●ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	● MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)

C11504	Solvent / Lubricant	● Isopropyle alcohol
C30618	Tube adaptor	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	● PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	● Raumedic-SIK 8806, Shore A 50-55
J1020091		
C10539	Protective cap	● PP (Polypropylene)
C11504	Solvent / Lubricant	● Isopropyle alcohol
C11511	Lubricant	● Aminofunctional Siloxane
C21791	Vented spike body	● PS (polystyrene) ● Colourant white, 50 % TiO ₂ in HIPS (high impact polystyrene)
J1040246		
C31536	Silicone fluid MED-420-350 cP	● Liquid copolymer
C31636	K-Zero Female Luer	● Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	● Liquid Silicone Rubber ● Phenyl Oil
Packaging		
	Tubular foil for packaging	● Film 500m x 350 mm x 40 µm ● Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Adhesive tape paper	● gummed paper tapes Polyeste filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	● HDPE (high density polyethylene)
	White paper for IFU printing	● 100% wood pulp
	Roll Label 105 x 148 mm	● Paper (blank white) and Adhesive
	Carton 400 x 265 x 205 mm	● Cardboard

	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Paper Volumat Line ON42 OP NF AirBlock	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	●100% Ethylene oxide, liquefied grade one (GAZ 4-134) ●alternative: > 99.9% Ethylene oxide