

Article code	M46444370
Article name acc. DoC	Infusion Set VL PN 95
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10478US	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C10845	filter (tape down clip = back of C10842)	●modified polyacryl
C10859	infusion flat filter 0.2 µm, neutrally charged	●Housing: polyacryl ●Filter: Supor® membrane, pore size 0.2 µm, hydrophilic polyethersulfone ●Vent membrane: 0.02 µm supported by PTFE (Polytetrafluoroethylene)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C30629	Tube adaptor	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) ●carrier material MABS (Terlux) with organic colour pigments (yellow)
C30911	Tube 2.94 x 4.02 mm	●PP (Polypropylene)
C30947	Tube 3 x 4.1 mm	●PUR (Polyurethane) + yellow UV stabilized Masterbatch PUR yellow
J1010116		
C10178	Drip chamber upper part 21 E	●PS (polystyrene) ●Colourant white, 50 % TiO ₂ in HIPS (high impact polystyrene)
C10221	moulded ring white for drip chamber 21	●PS (Polystyrene) ●Colourant, carrier material PS (Polystyrene), contains 50% TiO ₂
C10265	Drip chamber lower part 21 V92	●SBS (Styrene Butadiene Copolymer)
C10345	air ventilation flap	●1-propene homopolymer: PP (Polypropylene) ●1-propene copolymer (heterophasic copolymer)
C10539	Protective cap	●PP (Polypropylene)
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluoroethylene) Membrane ●Membrane Support Type: Non-woven Polyester
C11504	Solvent / Lubricant	●Isopropyle alcohol

Article code	Article name	Material
C11511	Lubricant	●Aminofunctional Siloxane
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch ●masterbatch orange
C10412	Roller clamp wheel	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch
J1010915		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30675	Revolving Luer Lock	●ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C30618	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	●PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	●Raumedic-SIK 8806, Shore A 50-55
J1040246		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●colour batch light blue
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent

Article code	Article name	Material
	Adhesive tape paper	●gummed paper tapes Polyeste filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	White paper for printing, IFU	●100% wood pulp
	Carton 400 x 265 x 205 mm	●Cardboard
	Carton board 390 x 200 mm	●Outliner / Fluting / Inliner
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Paper for VL PN95	●printed primary packaging paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide