

Article code	M46444600S
Article name acc. DoC	VL ON70 Oncology, PP material
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10478US	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C10605	connector Luer Lock	●MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene) - alternativ: C90106 MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)
C10845	filter (tape down clip = back of C10842)	●modified polyacryl
C10859	infusion flat filter 0.2 µm, neutrally charged	●Housing: polyacryl ●Filter: Supor® membrane, pore size 0.2 µm, hydrophilic polyethersulfone ●Vent membrane: 0.02 µm supported by PTFE (Polytetrafluoroethylene)
C10942	Tube 3.0 x 4.1 mm	●Polyether Polyurethane
C10942A	Polyurethane Tube Ø3,00 x Ø4,1 mm	●Polyether Polyurethane
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
	Solvent	●MEK (Methylethylketone)
C30911	Tube 2.94 x 4.02 mm	●PP (Polypropylene)
J1010110		
C10178	Drip chamber upper part 21 E	●PS (polystyrene) ●Colourant white, 50 % TiO₂ in HIPS (high impact polystyrene)
C10221	moulded ring white for drip chamber 21	●(98%) PS (Polystyrene) ●(2%) Colourant, carrier material PS (Polystyrene), contains 50% TiO₂
	Drip chamber lower part 21 V92	●SBS (Styrene Butadiene Copolymer)

Article code	Article name	Material
C10345	air ventilation flap	homopolymer: PP (Polypropylene) propene copolymer (heterophasic copolymer)
C10539	Protective cap	PP (Polypropylene)
C10808	Filter 15 µm drip chamber 21	●Filter: PA 6.6 (Polyamide) = Nylon Filter ●housing: PA 6 (Polyamide)
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ●Membrane Support Type: Non-woven Polyester
C11504	Solvent / Lubricant	●Isopropyle alcohol
	Lubricant	●Aminofunctional Siloxane
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch; ●masterbatch orange; Alternative to C90348 for FKN only:
C10412	Roller clamp wheel	●2.5% PS (Polystyrene) with silicone slip-masterbatch; Alternative to C90348 for FKN only: ●PS (Polystyrene)
J1010938		
C10477	Rotation clamp base	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene), alt. C90106
	Solvent / Lubricant	●Isopropyle alcohol
C30618	Tube adaptor	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butad
C30673	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
	Safety ring	●PP (Polypropylene)

Article code	Article name	Material
C30946	SIK-Tube Ø3,41xØ5,11	●Raumedic-SIK 8806, Shore A 50-55
Packaging		
C11003	PE bag	●LDPE (Low density polyethylene)
C11020	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
C11096	Edge protection	●corrugated cardboard
C11517	Adhesive tape paper	●gummed paper tapes Polyeste filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
C11945	White paper for IFU printing	●100% wood pulp
C12182	Roll Label 105 x 148 mm	●Paper (blank white) and Adhesive
C22614	perforatet film	●N/A
C31025A	Carton 400 x 265 x 205 mm	●Cardboard
C31045	Stretch film	●LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	●LDPE (low density polyethylene)
C31301	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
●THF (Tetrahydrofuran), also 1,4-Epoxy butane	Paper VP 133 x 267 mm	●Paper
C64571	Paper VP 500 Sets	
Sterilization		
C11520CN	Sterilization gas, Nanchang	●Ethylene oxide (90%) + Carbon dioxide (10%)
C11525	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
C31500	Sterilization gas, Poland	●100% Ethylene oxide, liquefied grade one (GAZ 4-134) ●alternative: > 99.9% Ethylene oxide