

Article code	M46444610
Article name acc. DoC	VL ON70 AirBlock PUR material, 0.2 micron filter
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C10478US	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C10605	connector Luer Lock	●MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene) - alternativ: MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)
C10845	filter	●modified polyacryl
C10859	infusion flat filter 0.2 µm, neutrally charged	●Housing: polyacryl ●Filter: Supor® membrane, pore size 0.2 µm, hydrophilic polyethersulfone ●Vent membrane: 0.02 µm supported by PTFE (Polytetrafluoroethylene)
C10942	Tube 3.0 x 4.1 mm	●Polyether Polyurethane
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C30883	Roller clamp safe	●Housing: PS (Polystyrene), colourant white ●Wheel: ABS (Acrylonitrile Butadiene Styrene), colourant blue
C30911	Tube 2.94 x 4.02 mm	●PP (Polypropylene)
J1010486		
C10221	moulded ring white for drip chamber 21	●(98%) PS (Polystyrene) ●(2%) Colourant, carrier material PS (Polystyrene), contains 50% TiO₂
C10345	air ventilation flap	●(50%) 1-propene homopolymer: PP (Polypropylene) ●(50%) 1-propene copolymer (heterophasic copolymer)
C10539	Protective cap	●PP (Polypropylene)
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ●Membrane Support Type: Non-woven Polyester

Article code	Article name	Material
C30879	Drip chamber upper part 21 safe	● PS (Polystyrene) ● (3%) carrier material S containing anorganic and organic colourants
C30881	Drip chamber lower part 21 safe	● SBS (Styrene Butadiene Copolymer)
C30882	Membrane safe 4 µm	● PES (Polyethersulfone)
J1010938		
C10477	Rotation clamp base	● MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	● Isopropyle alcohol
C30618	Tube adaptor	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	● PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	● Raumedic-SIK 8806
Packaging		
C11020	Tubular foil for packaging	● Film 500m x 350 mm x 40 µm ● Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
C11517	Adhesive tape paper	● gummed paper tapes ● Polyeste filament yarn for strengthening
C11521	PE-Band 90 x 30 x 0.05 mm	● HDPE (high density polyethylene)
C11945	White paper for IFU printing	● 100% wood pulp
C12182	Roll Label 105 x 148 mm	● Paper (blank white) and Adhesive
C31025A	Carton 400 x 265 x 205 mm	● Cardboard
C31045	Stretch film	● LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	● LDPE (low density polyethylene)
C31301	Compound foil 90µm	● Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
C64776	Paper Volumat Line ON70 AirBlock	● Paper
J1010486		
C11003	PE bag	● LDPE (Low density polyethylene)
Sterilization		
C11525	Indicator tape for EO sterilization	● Process indicator tape for Ethylene Oxide (EO) sterilization

Article code	Article name	Material
C31500	Sterilization gas	●100% Ethylene oxide, liquefied grade one