

Article code	M46444901Y
Article name acc. DoC	VL ON91 Light-protected, PUR material
Class acc. MDR 2017/745	IIa

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10478US	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C30629	Tube adaptor	● MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) ●carrier material MABS (Terlux) with organic colour pigments (yellow)
C30911	Tube 2.94 x 4.02 mm	●PP (Polypropylene)
C30947	Tube 3 x 4.1 mm	●PUR (Polyurethane) + 3,85% yellow UV stabilized Masterbatch PUR yellow
J1010111		
C10178	Drip chamber upper part 21 E	● PS (polystyrene) ● Colourant white, 50 % TiO ₂ in HIPS (high impact polystyrene)
C10221	moulded ring white for drip chamber 21	● (98%) PS (Polystyrene) ● (2%) Colourant, carrier material PS (Polystyrene), contains 50% TiO ₂
C10266	Drip chamber lower part 21 V92, yellow	● SBS (Styrene Butadiene Copolymer)(95%) ● masterbatch with UV-stabilizer, 5% yellow, Carrier: PS (Polystyrene) (5%)
C10345	air ventilation flap	●(50%) 1-propene homopolymer: PP (Polypropylene) ●(50%) 1-propene copolymer (heterophasic copolymer)
C10539	Protective cap	●PP (Polypropylene)
C10808	Filter 15 µm drip chamber 21	●Filter: PA 6.6 (Polyamide) = Nylon Filter ●housing: PA 6 (Polyamide)
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ●Membrane Support Type: Non-woven Polyester
C11504	Solvent / Lubricant	●Isopropyle alcohol
C11511	Lubricant	●Aminofunctional Siloxane
J1010810		

C10405	roller clamp housing, orange	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch; ●masterbatch orange ●PS (Polystyrene)
C10412	Roller clamp wheel	● PS (polystyrene); ●2.5% PS (Polystyrene) with silicone slip-masterbatch; ● PS (Polystyrene)
J1010915		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30675	Revolving Luer Lock	●ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C30618	Tube adaptor	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	●PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	●Raumedic-SIK 8806, Shore A 50-55
J1040246		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31636	K-Zero Female Luer	● Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		

	Tubular foil for packaging	● Film 500m x 350 mm x 40 µm ● Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	PE-Band 90 x 30 x 0.05 mm	● HDPE (high density polyethylene)
	White paper for IFU printing	● 100% wood pulp
	Roll Label 105 x 148 mm	● Paper (blank white) and Adhesive
	Adhesive tape	● BOPP (biaxially oriented polypropylene) film and rubber adhesive
	Carton 400 x 265 x 205 mm	● Cardboard
	Stretch film	● LLDPE (linear low density polyethylene)
	Stretch film foil plate	● LDPE (low density polyethylene)
	Compound foil 90µm	● Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Paper for M464444901Y, 133,67x267 mm	● Paper
Sterilization		
	Indicator tape for EO sterilization	● Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	● 100% Ethylene oxide, liquefied grade one (GAZ 4-134) ● alternative: > 99.9% Ethylene oxide