

Article code	M46445200S
Article name acc. DoC	Infusion Set VL PA 92
Class acc. 93/42/EEC	Ila

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10478	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C10478B	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C10607	LL neg. cone with tube fitment	●PVC (Polyvinylchloride) without plasticizer
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C11519	Solvent	●DCM (Dichlormethane) + 0.2 % Ethanole
C30710	NAC PRN needle free access device	●Housing: PC (Polycarbonate) ●Plug: Silicone Blue colourant - Pigmen mixture for silicone LubricanPoly(Methyl-3,3,3-Trifluoropropyl-siloxan)
C30727	Needle free access device	●Housing: PC (polycarbonate), Colourant and transparent part ● Stopper: Silicone
C30936	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) soft, noDOP, 34% TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30945	Tube 3.0 x 4.1 mm	●PVC (polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30945S	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C31567	Burette 150 ml with drip chamber	●Drip chamber: PVC (Polyvincylchloride) ●Top / Bottom cover burette: ABS (Acrylonitrile Butadiene Styrene) ●Cylinder tube: PETG (Polyethyleneterephthalate Glycol) ●Disk Filter: ABS (Acrylonitrile Butadiene Styrene) + PA (Polyamide) ●Air Filter: HDPE
C31580	solvent mixture	●61,3 % THF (Tetrahydrofuran) ●37,2 % MEK (Methylethylketone) ●1,5 % MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene); ●PS (Polystyrene) with silicone slip-masterbatch; masterbatch orange; ●PS (Polystyrene)

Article code	Article name	Material
C10412	Roller clamp wheel	●PS (polystyrene); ●2.5% PS (Polystyrene) with silicone slip-masterbatch; ●PS (Polystyrene)
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30605	Revolving Luer Lock	●luer cone: PVC (Polyvinylchloride) ●screwcap: ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C30618	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ● (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	●PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	●Raumedic-SIK 8806, Shore A 50-55
J1020091		
C10539	Protective cap	●PP (Polypropylene)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C11511	Lubricant	●Aminofunctional Siloxane
C21791	Vented spike body	●PS (polystyrene) ●Colourant white, 50 % TiO₂ in HIPS (high impact polystyrene)
Packaging		
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Adhesive tape paper	●gummed paper tapes - Polyeste filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)

Article code	Article name	Material
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	White paper for IFU printing	●100% wood pulp
	Carton 400 x 265 x 205 mm	●Cardboard
	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Paper for M46442200Y	●Paper
Sterilization		
	Sterilization gas, Nanchang	●Ethylene oxide (90%) + Carbon dioxide (10%)
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	●Ethylene oxide (100%) , liquefied grade one ● alt.: >99.9% Ethylene oxide