

Article code	M46445200Y
Article name acc. DoC	VL PA92 Paediatric, burette set
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10478	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C10478B	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C10607	LL neg. cone with tube fitment	●PVC (Polyvinylchloride) without plasticizer
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C11519	Solvent	●DCM (Dichlormethane) + 0.2 % Ethanol
C30936	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) soft, noDOP, 34% TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30945	Tube 3.0 x 4.1 mm	●PVC (polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30945S	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C31567	Burette 150 ml with drip chamber	●Drip chamber: PVC (Polyvinylchloride) ●Top / Bottom cover burette: ABS (Acrylonitrile Butadiene Styrene) ●Cylinder tube: PETG (Polyethyleneterephthalate Glycol) ●Disk Filter: ABS (Acrylonitrile Butadiene Styrene) + PA (Polyamide) ●Air Filter: HDPE
C31580	solvent mixture	●61,3 % THF (Tetrahydrofuran) ●37,2 % MEK (Methylethylketone) ●1,5 % MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene); ●PS (Polystyrene) with silicone slip-masterbatch; masterbatch orange; ●PS (Polystyrene)
C10412	Roller clamp wheel	●PS (polystyrene); ●2.5% PS (Polystyrene) with silicone slip-masterbatch; ●PS (Polystyrene)
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane

Article code	Article name	Material
C30605	Revolving Luer Lock	•luer cone: PVC (Polyvinylchloride) •screwcap: ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	•Isopropyle alcohol
C30618	Tube adaptor	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) • (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) • (3%) carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	•PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	•Raumedic-SIK 8806, Shore A 50-55
J1020091		
C10539	Protective cap	•PP (Polypropylene)
C11504	Solvent / Lubricant	•Isopropyle alcohol
C11511	Lubricant	•Aminofunctional Siloxane
C21791	Vented spike body	•PS (polystyrene) •Colourant white, 50 % TiO₂ in HIPS (high impact polystyrene)
J1040245		
C31536	Silicone fluid MED-420-350 cP	•Liquid copolymer
C31636	K-Zero Female Luer	•Eastman Tritan Copolyester MX731
C31637	K-Zero Male Luer	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) •colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	•Liquid Silicone Rubber •Phenyl Oil
J1040246		
C31536	Silicone fluid MED-420-350 cP	•Liquid copolymer
C31636	K-Zero Female Luer	•Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	•MABS (Methymethacrylate Acrylonitrile Butadiene Styrene) •colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	•Liquid Silicone Rubber •Phenyl Oil
J1040330		

Article code	Article name	Material
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane,
C11515	Solvent	●MEK (Methylethylketone)
C30605	Revolving Luer Lock	●luer cone: PVC (Polyvinylchloride) ●screwcap: ABS (Acrylonitrile Butadiene Styrene)
C30945	Tube 3.0 x 4.1 mm	●PVC (polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
C30945S	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31592	Tinopal Optical Brightener	●2,5-thiophenediylbis(5-tert-butyl-1,3-benzoxazole), color: Slight greenish yellow
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Adhesive tape paper	●gummed paper tapes - Polyeste filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	White paper for IFU printing	●100% wood pulp
	Carton 400 x 265 x 205 mm	●Cardboard
	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
	Paper for M46442200Y	●Paper
Sterilization		
	Sterilization gas, Nanchang	●Ethylene oxide (90%) + Carbon dioxide (10%)
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas, Poland	●Ethylene oxide (100%) , liquefied grade one ● alt.: >99.9% Ethylene oxide