

Article code	M46445475
Article name acc. DoC	VL SP 40-11
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10478US	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C10942	Tube 3.0 x 4.1 mm	●Polyether Polyurethane
C10942A	Polyurethane Tube Ø3,00 x Ø4,1 mm	●Polyether Polyurethane
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C30615	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30706	Back check valve	●Solvent: TEF (Tetrahydrofuran) ●Check Valve: MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene), Masterbatch ●Valve membrane: Silicone
C30911	Tube 2.94 x 4.02 mm	●PP (Polypropylene)
J1010118US		
C10178	Drip chamber upper part 21 E	●PS (polystyrene) ●Colourant white, 50 % TiO ₂ in HIPS (high impact polystyrene)
C10221	moulded ring white for drip chamber 21	●PS (Polystyrene) ●Colourant, carrier material PS (Polystyrene), contains 50% TiO ₂
C10265	Drip chamber lower part 21 V92	●SBS (Styrene Butadiene Copolymer)
C10345	air ventilation flap	●1-propene homopolymer: PP (Polypropylene) ●1-propene copolymer (heterophasic copolymer)
C10539	Protective cap	●PP (Polypropylene)
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ●Membrane Support Type: Non-woven Polyester
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch ●masterbatch orange

Article code	Article name	Material
C10412	Roller clamp wheel	●PS (polystyrene); ●PS (Polystyrene) with silicone slip-masterbatch
J1010915		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30675	Revolving Luer Lock	●ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C30618	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	●PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	●Silicone rubber, Shore A 50-55 ●without addition of colour pigment
J1040246		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●colour batch light blue
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		
C11020	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
C11517	Adhesive tape paper	●gummed paper tapes - Polyeste filament yarn for strengthening
C11521	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
C11600	Label sticker, size: 104 x 148 mm	●Paper for laserprinting

Article code	Article name	Material
C31025A	Carton 400 x 265 x 205 mm	●Cardboard, Flute C
C31043	Carton plate 390 x 200 mm	●Cardboard, Flute C
C31045	Stretch film	●LLDPE (linear low density polyethylene)
C31046	Stretch film foil plate	●LDPE (low density polyethylene)
C31301	Compound foil 90µm	●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)
C64556	Paper	●Paper
Sterilization		
C11525	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
C31500	Sterilization gas	●Ethylene oxide