

|                                       |                         |
|---------------------------------------|-------------------------|
| <b>Article code</b>                   | M46445475               |
| <b>Article name acc. DoC</b>          | Infusion Set VL SP40-11 |
| <b>Class acc. Directive 93/42/EEC</b> | IIa                     |

| Article code | Article name                           | Material                                                                                                                                                   |
|--------------|----------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| C10463       | Roberts clamp, white                   | ●ABS (Acrylonitrile Butadiene Styrene)                                                                                                                     |
| C10478US     | Lock plate                             | ●POM (Polyoxymethylene)<br>●carrier material Ionomer with anorganic and organic colour pigments (blue)                                                     |
| C10942       | Tube 3.0 x 4.1 mm                      | ●Polyether Polyurethane                                                                                                                                    |
| C10942A      | Polyurethane Tube Ø3,00 x Ø4,1 mm      | ●Polyether Polyurethane                                                                                                                                    |
| C11503       | Solvent                                | ●THF (Tetrahydrofuran), also 1,4-Epoxy butane                                                                                                              |
| C11515       | Solvent                                | ●MEK (Methylethylketone)                                                                                                                                   |
| C30615       | Tube adaptor                           | ●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)                                                                                                 |
| C30706       | Back check valve                       | ●Solvent: TEF (Tetrahydrofuran)<br>●Check Valve: MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene), Masterbatch<br>●Valve membrane: Silicone       |
| C30727       | needle free access device              | ●Housing: PC (Polycarbonate), Colourant and transparent part<br>●Stopper: Silicone                                                                         |
| C30911       | Tube 2.94 x 4.02 mm                    | ●PP (Polypropylene)                                                                                                                                        |
| J1010118US   |                                        |                                                                                                                                                            |
| C10178       | Drip chamber upper part 21 E           | ●PS (polystyrene)<br>●Colourant white, 50 % TiO <sub>2</sub> in HIPS (high impact polystyrene)                                                             |
| C10221       | moulded ring white for drip chamber 21 | ●PS (Polystyrene)<br>●Colourant, carrier material PS (Polystyrene), contains 50% TiO <sub>2</sub>                                                          |
| C10265       | Drip chamber lower part 21 V92         | ●SBS (Styrene Butadiene Copolymer)                                                                                                                         |
| C10345       | air ventilation flap                   | ●1-propene homopolymer: PP (Polypropylene)<br>●1-propene copolymer (heterophasic copolymer)                                                                |
| C10539       | Protective cap                         | ●PP (Polypropylene)                                                                                                                                        |
| C10836       | ventilation filter 3µm                 | ●Membrane Type: Hydrophobic<br>●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane<br>●Membrane Support Type: Non-woven Polyester             |
| J1010810     |                                        |                                                                                                                                                            |
| C10405       | roller clamp housing, orange           | ●PS (polystyrene)<br>●PS (Polystyrene) with silicone slip-masterbatch;<br>●masterbatch orange;<br>Alternative to C90348 for FKN only:<br>●PS (Polystyrene) |
| C10412       | Roller clamp wheel                     | ●PS (polystyrene);<br>●PS (Polystyrene) with silicone slip-masterbatch;                                                                                    |

| Article code  | Article name                              | Material                                                                                                                                                                                                            |
|---------------|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| J1010915      |                                           |                                                                                                                                                                                                                     |
| C10564        | Protective cap LL pos. with 1.2 µm filter | <ul style="list-style-type: none"> <li>●Housing: PP (Polypropylene) + Masterbatch</li> <li>●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane</li> </ul>            |
| C30675        | Revolving Luer Lock                       | <ul style="list-style-type: none"> <li>●ABS (Acrylonitrile Butadiene Styrene)</li> </ul>                                                                                                                            |
| J1010938      |                                           |                                                                                                                                                                                                                     |
| C10477        | Rotation clamp base                       | <ul style="list-style-type: none"> <li>●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)</li> </ul>                                                                                                        |
| C11504        | Solvent / Lubricant                       | <ul style="list-style-type: none"> <li>●Isopropyle alcohol</li> </ul>                                                                                                                                               |
| C30618        | Tube adaptor                              | <ul style="list-style-type: none"> <li>●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)</li> <li>●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)</li> </ul> |
| C30638        | Tube adaptor                              | <ul style="list-style-type: none"> <li>●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)</li> </ul>                                                                                                        |
| C30673        | Tube adaptor                              | <ul style="list-style-type: none"> <li>●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)</li> <li>●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)</li> </ul> |
| C30674        | Safety ring                               | <ul style="list-style-type: none"> <li>●PP (Polypropylene)</li> </ul>                                                                                                                                               |
| C30946        | SIK-Tube Ø3,41xØ5,11                      | <ul style="list-style-type: none"> <li>●Raumedic-SIK 8806, Shore A 50-55</li> </ul>                                                                                                                                 |
| Packaging     |                                           |                                                                                                                                                                                                                     |
|               | Tubular foil for packaging                | <ul style="list-style-type: none"> <li>●Film 500m x 350 mm x 40 µm</li> <li>●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent</li> </ul>                          |
|               | Adhesive tape paper                       | <ul style="list-style-type: none"> <li>●gummed paper tapes</li> <li>●Polyeste filament yarn for strengthening</li> </ul>                                                                                            |
|               | PE-Band 90 x 30 x 0.05 mm                 | <ul style="list-style-type: none"> <li>●HDPE (high density polyethylene)</li> </ul>                                                                                                                                 |
|               | Label sticker, size: 104 x 148 mm         | <ul style="list-style-type: none"> <li>●Paper for laserprinting</li> </ul>                                                                                                                                          |
|               | Carton 400 x 265 x 205 mm                 | <ul style="list-style-type: none"> <li>●Cardboard</li> </ul>                                                                                                                                                        |
|               | Carton board 390 x 200 mm                 | <ul style="list-style-type: none"> <li>●Outliner / Fluting / Inliner</li> </ul>                                                                                                                                     |
|               | Stretch film                              | <ul style="list-style-type: none"> <li>●LLDPE (linear low density polyethylene)</li> </ul>                                                                                                                          |
|               | Stretch film foil plate                   | <ul style="list-style-type: none"> <li>●LDPE (low density polyethylene)</li> </ul>                                                                                                                                  |
|               | Compound foil 90µm                        | <ul style="list-style-type: none"> <li>●Multilayer flat film PA (Polyamide) /adhesive layer / PE (Polyethylene)</li> </ul>                                                                                          |
|               | Paper                                     | <ul style="list-style-type: none"> <li>●Paper</li> </ul>                                                                                                                                                            |
| Sterilization |                                           |                                                                                                                                                                                                                     |
|               | Sterilization gas                         | <ul style="list-style-type: none"> <li>●Ethylene oxide</li> </ul>                                                                                                                                                   |
|               | Indicator tape for EO sterilization       | <ul style="list-style-type: none"> <li>●Process indicator tape for Ethylene Oxide (EO) sterilization</li> </ul>                                                                                                     |