

Article code	M46445865
Article name acc. DoC	VL PR72-11, 0.2 micron filter
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10478	Lock plate	●POM (Polyoxymethylene) ●carrier material Ionomer with anorganic and organic colour pigments (blue)
C10845	filter (tape down clip)	●modified polyacryl
C10859	infusion flat filter 0.2 µm, neutrally charged	●Housing: polyacryl ●Filter: Supor® membrane, pore size 0.2 µm, hydrophilic polyethersulfone ●Vent membrane: 0.02 µm supported by PTFE (Polytetrafluoroethylene)
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
C30706	Back check valve	●Solvent: TEF (Tetrahydrofuran) ●Check Valve: MABS (Methymethacrylate Acrylonitrile Butadiene Styrene), Masterbatch ●Valve membrane: Silicone
C30945S	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) with TOTM (Trioctyl Trimellitate) plasticizer
C31580	solvent mixture	●THF (Tetrahydrofuran) ●MEK (Methylethylketone) ●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
J1010116		
C10178	Drip chamber upper part 21 E	●PS (polystyrene) ●Colourant white, 50 % TiO ₂ in HIPS (high impact polystyrene)
C10221	moulded ring white for drip chamber 21	●PS (Polystyrene) ●Colourant, carrier material PS (Polystyrene), contains 50% TiO ₂
C10265	Drip chamber lower part 21 V92	●SBS (Styrene Butadiene Copolymer)
C10345	air ventilation flap	●1-propene homopolymer: PP (Polypropylene) ●1-propene copolymer (heterophasic copolymer)
C10539	Protective cap	●PP (Polypropylene)

Article code	Article name	Material
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ●Membrane Support Type: Non-woven Polyester
C11504	Solvent / Lubricant	●Isopropyle alcohol
C11511	Lubricant	●Aminofunctional Siloxane
J1010831		
C10412	Roller clamp wheel	●PS (polystyrene); ●PS (Polystyrene) with silicone slip-masterbatch
C10432	Roller clamp housing	●PS (Polystyrene) ●PS (Polystyrene) with silicone ●Carrier material PS (Polystyrene) with anorganic and organic colour pigments (blue)
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30605	Revolving Luer Lock	●luer cone: PVC (Polyvinylchloride) ●screwcap: ABS (Acrylonitrile Butadiene Styrene)
J1010938		
C10477	Rotation clamp base	●MABS (Methymethacrylate Acrylonitrile Butadiene Styrene)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C30618	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30638	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30673	Tube adaptor	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●carrier material PS (Polystyrene) with anorganic and organic pigments and dyes (light green)
C30674	Safety ring	●PP (Polypropylene)
C30946	SIK-Tube Ø3,41xØ5,11	●Raumedic-SIK 8806, Shore A 50-55
J1040246		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31638	K-Zero Y-Site Body	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●colour batch light blue

Article code	Article name	Material
C31639	K-Zero Microbial Barrier	• Liquid Silicone Rubber • Phenyl Oil
Packaging		
	PE bag	• LDPE (Low density polyethylene)
	Tubular foil for packaging	• Film 500m x 350 mm x 40 µm • Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Adhesive tape paper	• gummed paper tapes • Polyester filament yarn for strengthening
	PE-Band 90 x 30 x 0.05 mm	• HDPE (high density polyethylene)
	Label sticker, size: 104 x 148 mm	• Paper for laserprinting
	PE bag with bottomseal, transparent	• LDPE (low density polyethylene) • Film thickness: 45 µm
	Carton 400 x 265 x 205 mm	• Cardboard
	Carton board 390 x 200 mm	• Outliner / Fluting / Inliner
	Stretch film	• LLDPE (linear low density polyethylene)
	Stretch film foil plate	• LDPE (low density polyethylene)
	PE (Polyethylene) film broadband	• PE (Polyethylene), white
	Compound foil 90µm	• Multilayer flat film PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Paper for VL PR 72-11	• Paper
Sterilization		
	Indicator tape for EO sterilization	• Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	• Ethylene oxide