

Article code	M77460026
Article name acc. DoC	Secondary Line FX PVC-free, with 0,2 µ Filter
Class acc. Directive 93/42/EEC	Is

Article code	Article name	Material
C10463	Roberts clamp, white	●ABS (Acrylonitrile Butadiene Styrene)
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C10605	connector Luer Lock	●MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)
C10845	filter (tape down clip)	●modified polyacryl
C10859	infusion flat filter 0.2 µm, neutrally charged	●Housing: polyacryl ●Filter: Supor® membrane, pore size 0.2 µm, hydrophilic polyethersulfone ●Vent membrane: 0.02 µm supported by PTFE (Polytetrafluoroethylene)
C10942 C10942A	Tube 3.0 x 4.1 mm	●Polyether Polyurethane
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C11515	Solvent	●MEK (Methylethylketone)
J1020017		
C10539	Protective cap	●PP (Polypropylene)
C21795	Spike	●PS (Polystyrene) ●carrier material PS (Polystyrene) containing anorganic and organic colourants
C21796	Freeflex claw	●ABS (Acrylonitrile Butadiene Styrene) ●PS (Polystyrene) with Inorganic and organic pigments
Packaging		
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Carton box	●Cardboard
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Adhesive tape paper	●gummed paper tapes Polyeste filament yarn for strengthening
	PE banderole 60 x25 x 0.10 mm	●LDPE (low density polyethylene)
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	White paper for IFU printing	●100% wood pulp

Article code	Article name	Material
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Paper Secondary line 100 x 180 mm	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide