

Article code	M78401045Y
Article name acc. DoC	K-Zero Extension 10cm
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C10425	Slide clamp, blue	●ABS (Acrylonitrile Butadiene Styrene) ●Carrier: ABS (Acrylonitrile Butadiene Styrene) - Components: Inorganic and organic pigments
C10603	male connector Luer Lock	●MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene) ●MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)
C10636	Female connector LL	●PVC (Polyvinylchloride) without plasticizer
C11519	Solvent	●DCM (Dichlormethane) + 0.2 % Ethanole
C11553	Solvent	●CHN (Cyclohexanone)
C30937	No-DOP-Schlauch 1,5 x 2,7 mm	●PVC (Polyvinylchloride) containing plasticizer
C31640	K-Zero Vented Male Luer End Cap	●LDPE (low density polyethylene)
J1040245		
C31536	Silicone fluid MED-420-350 cP	●Liquid copolymer
C31636	K-Zero Female Luer	●Eastman Tritan Copolyester MX731
C31637	K-Zero Male Luer	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene) ●colour batch light blue, 3%
C31639	K-Zero Microbial Barrier	●Liquid Silicone Rubber ●Phenyl Oil
Packaging		
	Carton 400 x 220 x 380 mm	●Corrugated cardboard, Flute B
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Adhesive tape paper	●gummed paper tapes Polyeste filament yarn for strengthening
	Label sticker, size: 104 x 148 mm	●Paper for laser printing
	White paper for printing, IFU	●100% wood pulp
	Stretch film	●LLDPE (linear low-density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Paper for M78401045Y	●Paper
Sterilization		

Article code	Article name	Material
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide