

Article code	M93000030
Article name acc. DoC	Injectomat® Syringe 50 ml, light protected
Class acc. MDR 2017/745	Ila

Article code	Article name	Material
C11511	Lubricant	●Aminofunctional Siloxane
C31545	Print color black	●Ink
C31644	Live 360 Medical Fluid, 12,500 cSt	●silicone fluid
C31645	Liveo Q7-9180 Silicone Fluid, 1 cSt	●silicone fluid
C31725	Stopper for syringe	●Polyisoprene rubber + 2% Carbon Black
C31726	Plunger 50ml	●PP (Polypropylene) transparent
C31728	Syringe barrel 50ml orange	●PP (Polypropylene) ●colourant orange
Packaging		
	Roll Label 105 x 148 mm	●Paper (blank white) and Adhesive
	Adhesive tape	●BOPP (biaxially oriented polypropylene) film and rubber adhesive
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Carton for syringes 50ml	●Cardboard
	Paper	●Gas Paper 70 HP (Medical paper)
	Fitform Foil 150 µm	●PA/PE (Polyamide) / Polyethylene
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide