

Article code	F2859311
Article name acc. DoC	Freka® Mix Bag 3000 mL + 3-filling leads
Class acc. MDR 2017/745	Is

Article code	Article name	Material
311-01	Tube foil 0.35 mm	●EVA (ethylene vinyl acetate copolymer)
322041	Coex tube 6.3 x 7.95 x 40 mm	●EVA (ethylene vinyl acetate copolymer)/PUR (Polyurethane)
322042	Coex tube 4.9 x 6.5 mm	●EVA (ethylene vinyl acetate copolymer)/PUR (Polyurethane)
330152	Breakable injection site with protective cap	●ABS (Acrylonitrile Butadiene Styrene) ●PP (Polypropylene) ●Polyisoprene (Latex-free)
330177	Non-reopening clamp, nature	●PP (Polypropylene)
330251	Locking cap female, blue	●ABS(Acrylonitrile-Butadiene-Styrene) ●Colourant
330254	LL connector male, white	●ABS (Acrylonitrile Butadiene Styrene)
330255	Infusion closure (latex-free, aluminium-free)	●PC (Polycarbonate), Chlorobutyle
330289	Handle board, blue	●PP (Polypropylene), Colourant
330305 + 330306	Injection port for Breakable Injection Site + Protective cap for injection port	●PC (Polycarbonate)/PP (Polypropylene), Polyisoprene rubber ●LDPE (Low-density Polyethylene)
350-01	Bag imprint	●Stamp foil blue
Material 239	Solvent	●THF (Tetrahydrofuran)
J1040332		
C11553	Solvent	●CHN (Cyclohexanone)
C31740	Luer Lock connector female, white	●ABS (Acrylonitrile Butadiene Styrene)
C31747	Clamp, white	●PP (Polypropylene)
C31749	4-way Y-Connector	●PVC (Polyvinylchloride), plasticizer TEHTM (Tri-2-(Ethylhexyl) Trimellitate)
C31754	Spike Assembly, Filter 3µm, Cap	●Spike: ABS (Acrylonitrile Butadiene Styrene); ●Filter 3 µm: PE (Polyethylene) / Versapor 3000RC; ●Cap: LDPE (Low-Density Polyethylene)
C60963	Tube 4.8 x 6.8	●PVC (Polyvinylchloride) with plasticizer TEHTM (Tri-(2-ethylhexyl) trimellitate)
Packaging		
340100	Sterile gas paper	●Paper
340102	Compound foil 80 µm	●PA (Polyamide) / PE (Polyethylene)

Article code	Article name	Material
340104	Foil band 180 x 100 mm	●HDPE (high density polyethylene)
341020	Carton OD 400 x 400 x 250 mm	●corrugated cardboard
342014	Label sticker 105 x 148 mm	●Paper
350-04	Transfer film	●Transfer film, black
350-05	Transfer print	●Transfer foil, black
350105	Flexoprint	●Printing ink, black
Sterilization		
322036	Sterilization indicator, brown	●Sterilization indicator ETO
Material 240	Sterilization	●ETO (Ethylene Oxide)