

Material List – Vasodrop Safety

Article code	See attachment I
Article name acc. DoC	Vasodrop Safety
Class acc. Directive 93/42/EEC	IIa

Article name	Material
Luer Lock Male Cap	Polypropylene
Vent Fitting	ABS
Filter	High density polyethylene
Needle hub	ABS
Catheter hub with wings and port	Polypropylene
Valve	Silicone tube
Injection site protector	Polypropylene
Caulking pin	Polyoxymethylene
Tube catheter	PUR + BaSO ₄
Needle	Stainless steel
Needle protector	Low density polyethylene
Safety housing	Polyamide
Safety fingers	Polysulfone
Sliding bush	Polycarbonate
Needle blocker	Stainless steel
Packaging	
Primary	Thermoformed and thermosealed blister pack, with medical grade mesh paper and plastic film
Secondary	Cellulose cardboard box
For transport	Corrugated cardboard
Sterilization	
Sterilization gas	ETO (Ethylene Oxide)

Annex 1 Covered products:

Article code	Article name
8301700	Vasodrop Safety 2 14G 45mm
8301701	Vasodrop Safety 2 16G 45mm
8301702	Vasodrop Safety 2 17G 45mm
8301704	Vasodrop Safety 2 18G 45mm
8301703	Vasodrop Safety 2 18G 32mm
8301705	Vasodrop Safety 2 20G 32mm
8301706	Vasodrop Safety 2 22G 25mm
8301707	Vasodrop Safety 2 24G 19mm