

Article code	Z072025F (M78402060)
Article name acc. DoC	Extension Y BCV
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C10479	Roberts Clamp blue	● PP (Polypropylene) ● colour batch: PE (Polyethylene)
C10511	protective cap LL	● PS (polystyrene) with butadien HIPS (high impact polystyrene) ● PS (polystyrene) with silicone
C10557	Protective cap, transparent	● PP (Polypropylene)
C10644	Kone LL neg. for Tube 1,0 x 2,0	● PVC (Polyvinylchloride) without plasticizer
C10664	Kone LL pos. f. Innenver-	● PVC (Polyvinylchloride) without plasticizer
C10770	Back check valve	● Housing upper part: SAN (Styrene acrylonitrile copolymer) ● Housing lower part: SAN (Styrene acrylonitrile copolymer) ● Colourant: white ● Valve: Silicone rubber
C11519	Solvent	● DCM (Dichlormethane) + 0.2 % Ethanol
C11553	Solvent	● CHN (Cyclohexanone)
C30721	Y-connector	● MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30900	PVC-Schlauch 1,0 x 2,0 mm	● PVC (Polyvinylchloride) with TOTM
C30936	Tube 3.0 x 4.1 mm	● PVC (Polyvinylchloride) soft, noDOP, 34% TEHTM (Tri-(2-ethylhexyl) trimellitate)
J1010920		
C10505	Protective cap lock	● PP (Polypropylene)
C10605	connector Luer Lock	● MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)
Packaging		
	Compound foil 70 µm	● PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Adhesive tape paper	● gummed paper tapes ● Polyester filament yarn for strengthening
	PE banderole 60 x25 x 0.10 mm	● LDPE (low density polyethylene)
	Label sticker, size: 97 x 67.7 mm	● Paper 70 g / qm ● adhesive (Acrylate-dispersion)
	White paper for IFU printing	● 100% wood pulp
	Carton box 195 x 195 x 110 mm	● Corrugated cardboard, color brown inside and outside
	PE bag with bottomseal, transparent	● LDPE (low density polyethylene) ● Film thickness: 45 µm
	Stretch film	● LLDPE (linear low density polyethylene)
	Stretch film foil plate	● LDPE (low density polyethylene)

Article code	Article name	Material
	PE (Polyethylene) film broadband	●PE (Polyethylene), white
	Paper "Extension"	●Paper
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide