

Article code	Z072039F (M78402070)
Article name acc. DoC	Extension ASV Y BCV
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C10479	Roberts Clamp blue	●PP (Polypropylene) ●colour batch: PE (Polyethylene)
C10511	protective cap LL	●PS (polystyrene) with butadien HIPS (high impact polystyrene) ●PS (polystyrene) with silicone
C10557	Protective cap, transparent	●PP (Polypropylene)
C10644	Kone LL neg. for Tube 1,0 x 2,0	●PVC (Polyvinylchloride) without plasticizer
C10664	Kone LL pos. f. Innenver-	●PVC (Polyvinylchloride) without plasticizer
C10770	Back check valve	●Housing upper part: SAN (Styrene acrylonitrile copolymer) ●Housing lower part: SAN (Styrene acrylonitrile copolymer) ●Colourant: white ●Valve: Silicone rubber
C10786	Anti-Syphon Valve	●Upper Part: PC (Polycarbonate) ●Lower Part: PC (Polycarbonate) ●Internal Part: Silicone
C11519	Solvent	●DCM (Dichlormethane) + 0.2 % Ethanol
C11553	Solvent	●CHN (Cyclohexanone)
C30721	Y-connector	●MABS (Methylmethacrylate Acrylonitrile Butadiene Styrene)
C30900	PVC-Schlauch 1,0 x 2,0 mm	●PVC (Polyvinylchloride) with TOTM
C30936	Tube 3.0 x 4.1 mm	●PVC (Polyvinylchloride) soft, noDOP, 34% TEHTM (Tri-(2-ethylhexyl) trimellitate)
J1010920		
C10505	Protective cap lock	●PP (Polypropylene)
C10605	connector Luer Lock	●MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)
Packaging		
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	Adhesive tape paper	●gummed paper tapes ●Polyeste filament yarn for strengthening
	PE banderole 60 x25 x 0.10 mm	●LDPE (low density polyethylene)
	Label sticker, size: 97 x 67.7 mm	●Paper 70 g / qm ●adhesive (Acrylate-dispersion)
	White paper for IFU printing	●100% wood pulp
	Carton box 195 x 195 x 110 mm	●Corrugated cardboard, color brown inside and outside

Article code	Article name	Material
	PE bag with bottom seal, transparent	●LDPE (low density polyethylene) ●Film thickness: 45 µm
	Stretch film	●LLDPE (linear low-density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	PE (Polyethylene) film broadband	●PE (Polyethylene), white
	Paper "Extension"	●Paper
Sterilization		
	Sterilization gas, Nanchang	●Ethylene oxide
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization