

Article code	Z072710FCB
Article name acc. DoC	VS 10
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C30910	Tube 2.94 x 4.09 mm	●PVC (Polyvinylchloride), with plasticizer DEHP
J1010412		
C10110	Drip chamber upper part	●PS (Polystyrene)
C10209	Drip chamber lower part 18	●SBS (Styrene Butadiene Copolymer)
C10216	Moulded ring white for drip chamber 18	●PS (Polystyrene) ●Colourant, carrier material PS (Polystyrene), contains 50% TiO2
C10345	air ventilation flap	●1-propene homopolymer: PP (Polypropylene) ●1-propene copolymer (heterophasic copolymer)
C10539	Protective cap	●PP (Polypropylene)
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ●Membrane Support Type: Non-woven Polyester
C10838	Disc filter 15 µm for US	●PA 6.6 (Polyamid 6.6)
C11504	Solvent / Lubricant	●Isopropyle alcohol
C11511	Lubricant	●Aminofunctional Siloxane
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch; ●masterbatch orange
C10412	Roller clamp wheel	●PS (polystyrene); ●PS (Polystyrene) with silicone slip-masterbatch
J1010920		
C10505	Protective cap lock	●PP (Polypropylene)
C10605	connector Luer Lock	●MABS (Methylmetacrylate Acrylonitrile Butadiene Styrene)
Packaging		
	Tubular foil for packaging	●Film 500m x 350 mm x 40 µm ●Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Carton 400 x 265 x 380 mm	●Corrugated cardboard
	Compound foil 70 µm	●PA (Polyamide) / adhesive layer / PE (Polyethylene)
	PE-Band 90 x 30 x 0.05 mm	●HDPE (high density polyethylene)
	Label sticker, size: 104 x 148 mm	●Paper for laserprinting
	White paper for IFU printing	●100% wood pulp

Article code	Article name	Material
	Adhesive tape	●BOPP (biaxially oriented polypropylene) film and rubber adhesive
	Stretch film	●LLDPE (linear low density polyethylene)
	Stretch film foil plate	●LDPE (low density polyethylene)
	Paper VS 10	●Paper
	Carton 1185 x 785	●Carton
	PE-bag with bottemseal, transparent	●LDPE (low density polyethylene) ●Film thickness: 80 µm
Sterilization		
	Indicator tape for EO sterilization	●Process indicator tape for Ethylene Oxide (EO) sterilization
	Sterilization gas	●Ethylene oxide